

# **Progenesis Q1 for proteomics User Guide**

**Analysis workflow guidelines for DDA data**

# **Waters**

**THE SCIENCE OF WHAT'S POSSIBLE.™**

## Contents

Introduction.....	3
How to use this document .....	3
How can I analyse my own runs using Progenesis QI for proteomics? .....	3
LC-MS Data used in this user guide .....	3
Workflow approach to LC-MS run analysis.....	4
Restoring the Tutorial .....	5
Stage 1: Import Data and QC review of LC-MS data set.....	6
Stage 2A: Automatic Processing of your data .....	7
Stage 2B: After Automatic Processing.....	12
Stage 3: Licensing.....	13
Stage 4: Review Alignment .....	14
Reviewing Quality of Alignment.....	15
Stage 5A: Filtering.....	18
Stage 5B: Reviewing Normalisation .....	22
Stage 6: Experiment Design Setup for Analysed Runs .....	25
Stage 7: Review Peak Picking and editing of results.....	28
Stage 8: Peptide Ion Statistics on selected peptide ions .....	37
Stage 9: Identify peptides .....	41
Stage 10: QC Metrics .....	45
Stage 11: Refine Identifications.....	47
Resolving Conflicts.....	49
Protein Grouping .....	49
Protein Quantitation options .....	51
Stage 12: Review Proteins .....	54
Stage 13: Exporting Protein Data .....	61
Exporting Protein Data to Pathways Tool(s).....	62
Stage 14: Exporting identifications for submission to PRIDE .....	65
Stage 15: Protein Statistics .....	67
Stage 16: Reporting.....	68
Creating an Inclusion list .....	70
Appendix 1: Stage 1 Data Import and QC review of LC-MS data set.....	72
Appendix 2: Stage 1 Processing failures .....	76
Appendix 3: Stage 1 Data QC review and addition of exclusion areas.....	78
Appendix 4: Licensing runs (Stage 3).....	79
Appendix 5: Manual assistance of Alignment .....	80
Appendix 6: Within-subject Design.....	84
Appendix 7: Power Analysis (Peptide Ion Stats).....	86
Appendix 8: Resolve Conflicts .....	87
Appendix 9: Using Clip Gallery to Save and Export Pictures and Data .....	91
Appendix 10: Exporting protein identities to IPA.....	93
Appendix 11: Exporting protein identities to MetaCore.....	96
Appendix 12A Search engine parameters for Mascot .....	98
Appendix 12B Use Additional Export Fields in Mascot for PRIDE .....	99
References .....	99

## Introduction

This user guide takes you through a complete analysis of 6 LC-MS runs with 2 groups (3 replicate runs per group) using the unique Progenesis QI for proteomics workflow. It starts with LC-MS data file loading then Alignment, followed by Analysis that creates a list of interesting peptide ions (peptides) which are explored within Peptide Ion Stats using multivariate statistical methods then onto Protein identity and Reporting.

To allow ease of use the tutorial is designed to start with the restoration of an Archived experiment where the data files have already been loaded. The document covers all the stages in the workflow, initially focusing on the Automatic Processing of the data then the use of the individual analysis stages.

If you are using your own data files please refer to Appendix 1 (page 72) then start at page 7.

## How to use this document

You can print this user guide to help you work hands-on with the software. The complete user guide takes about 60 to 90 minutes (dependant on PC spec) and is divided into two sections. This means you can perform the first half focused on LC-MS run alignment and analysis then complete the second half of analysis exploring comparative differences and Protein identity at a convenient time.

If you experience any problems or require assistance, please contact us at [support@nonlinear.com](mailto:support@nonlinear.com)

## How can I analyse my own runs using Progenesis QI for proteomics?


You can freely explore the quality of your LC-MS data using Data Import and then licence your own LC-MS runs using this evaluation copy of Progenesis QI for proteomics. Instructions on how to do this are included in a section at the end of the user guide document. Alternatively if you would like to arrange a demonstration in your own laboratory contact [support@nonlinear.com](mailto:support@nonlinear.com) and we will help you.

## LC-MS Data used in this user guide

NLD would like to thank Dr Robert Parker and Prof Haroun Shah at the Health Protection Agency, London, UK for providing the example data used in this user guide as well as invaluable discussion on the handling of the data.

## Workflow approach to LC-MS run analysis

Progenesis QI for proteomics adopts an intuitive **Workflow** approach to performing comparative LC-MS data analysis. The following user guide describes the various stages of this workflow (see below) focusing mainly on the stages from Alignment to Report.




Stage	Description	Page
	<b>LC-MS Import Data:</b> Selection and review of data files for analysis	6
	<b>Automatic Processing:</b> setting up steps for automatic processing	7
	<b>After Automatic Processing:</b> how to work with auto analysed data	12
	<b>Licensing:</b> allows licensing of individual data files when there is no dongle attached (Appendix 4 (page 79))	13
	<b>Review Alignment:</b> automatic and manual LC-MS run alignment	14
	<b>Filtering:</b> defining filters for peaks based on Retention Time, m/z, Charge State and Number of Isotopes.	18
	<b>Review Normalisation:</b> exploring LC-MS normalisation	22
	<b>Experiment Design Setup:</b> defining one or more group set ups for analysed aligned runs	25
	<b>Review Peak Picking:</b> review and validate results, edit peak detection, tag groups of peaks and select peaks for further analysis	28
	<b>Peptide Ion Statistics:</b> performing multivariate statistical analysis on tagged and selected groups of peptide ions	37
	<b>Identify Peptides:</b> managing export of MS/MS spectra to, and import of peptide ids from Peptide Search engines.	41
	<b>QC metrics:</b> quality control charts for experimental/analysed data	45
	<b>Refine Identifications:</b> manage filtering of peptide ids	47
	<b>Resolving Conflicts:</b> validation and resolution of peptide id conflicts for data entered from Database Search engines	49
	<b>Review proteins:</b> review protein and peptide identity and data export	54
	<b>Protein Statistics:</b> multivariate statistical analysis on proteins	67
	<b>Report:</b> generate a report for proteins and/or peptides	68

## Restoring the Tutorial


Open Progenesis QI for proteomics and download the Compressed (.zip) Tutorial Archive file from the 'User guide and tutorial' link shown below, placing it in a **new folder** on your desktop. Before restoring the tutorial in the software **you must** first right click on the (.zip) file and extract it to the same folder.

Now you can restore the uncompressed Tutorial archive file. To do this, first locate the **Progenesis QI.P\_Tutorial DDA.Progenesis QIP Archive** file using the **Open** button and press Open.




This opens the 'Import Experiment from archive' dialog.

Select the **Create a new experiment** option and select the folder in which you placed the archive, using Browse.



Then press **Import**.




**Note:** use the **Replace an existing experiment** option if you want to over-write an existing version of the tutorial.

**Tip:** at each stage in the software there are links to more information and help on the website.


## Stage 1: Import Data and QC review of LC-MS data set

The tutorial will now open at the **Import Data** stage (see below).




Each data file appears as a 2D representation of the run.

**Note:** the **Experiment Properties** are available from the File menu. These were selected when the experiment was created (see Appendix 1, page 72). The tutorial data is profile data.



**Tip:** the **'Mask areas for peak picking'** facility allows you to examine and exclude areas (usually early and/or late in the LC dimension (Retention Time)) that appear excessively noisy due to capture of data during column regeneration (see Appendix 3, page 78). This is not required for this data set.



**Note:** use the **Remove Run** to remove run(s) from the current experiment.

Now start the Automatic Processing.

## Stage 2A: Automatic Processing of your data

The Automatic Processing of your data can be set up and started before the import of your data has been completed by clicking on **Start automatic processing**.

**Note:** for this tutorial the data has been imported already.

**2 Perform automatic processing**  
Click the button below to start automatic processing. This will analyse your data as far as possible, before suggesting the next step.

**Start automatic processing**

Automatic processing can be started while runs are importing.


[Learn more about automatic processing](#)

Setup of processing steps in the Analysis Workflow, up to and including Identify Peptides, can be performed in the Automatic Processing wizard, these include:

- Automatic alignment of all runs to a reference run
- Automatic peak picking
- Creating an Experiment design
- Choosing a Quantitation method
- Identification of peptides (only available for MSe and HDMSe data formats)

In this tutorial example you have 6 Data Dependently Acquired (DDA) LC-MS runs, so the automatic identification of peptides is unavailable in the automatic processing.

As the runs have already been imported, click **Start automatic processing** to setup the Automatic processing wizard.



**Progenesis QI for proteomics provides three methods** for choosing the alignment reference run:

### 1. Assess all runs in the experiment for suitability

This method compares every run in your experiment to every other run for similarity, then selects the run with the greatest similarity to all other runs as the alignment reference.

If you have no prior knowledge about which of your runs would make a good reference, then this choice will normally produce a good alignment reference for you. This method, however, can take a long time for a large number of runs.

### 2. Use the most suitable run from candidates that I select

This method asks you to make a selection of reference candidates; the automatic algorithm then chooses the best reference from this subset of runs.

This method is appropriate when you have some prior knowledge of your runs suitability as references: i.e when all the candidate runs are pooled samples or when all the candidates are from a condition that displays the largest set of common peptide ions.

### 3. Use this run


This method allows you to manually choose the reference run.

Manual selection gives you full control, but there are a couple of risks to note:

- If you choose a pending run which subsequently fails to load, alignment will not be performed.
- If you choose a run before it fully loads, and it turns out to have chromatography issues, alignment will be negatively affected (for this reason we recommend that you let your reference run fully load and assess its chromatography before loading further runs).


**For this tutorial** we will select the first option, (See Appendix 1, page 72 for more details on using the other options).

You will now be asked if you want to align your runs automatically.



The default is for automatic alignment, click **Next**.



The next page of the processing wizard asks you if you want to **Perform peak picking** and allows you to set appropriate parameters.



For the purposes of this User guide we will use the default settings for peak picking except in the case of Retention Time limits where we will limit the peak picking to between 10 and 75 min.




**Note:** for more details on setting Peak Picking parameters refer to the section on Filtering (page 18) Click on **Set parameters** and select the fourth tab to set **Retention time limits** for the detection. The default limits are displayed in the unticked boxes and detection will use these values whenever the box(s) are unticked. Enter values of 10 and 75 min and tick the boxes as shown below.



Click **OK** to return to the Automatic Processing Wizard and click **Next** to Define an Experiment design.

To handle the grouping of your run data you can make use of **sample tracking information** that has been stored in a spreadsheet at the time of sample collection and/or preparation.




**Note:** if a spreadsheet file of your data is not available you can create your experiment designs after the automatic analysis of the runs

For this example there is a **Tutorial Groups.csv** file available with the Tutorial Archive you restored at the beginning of this tutorial exercise.

Sample Name	Conditions	Date of Collection	Location
A1	A	02/03/2011	Fridge A
A2	A	02/03/2011	Fridge A
A3	A	02/03/2011	Fridge B
C1	C	06/03/2011	Fridge B
C2	C	06/03/2011	Fridge B
C3	C	06/03/2011	Fridge A

Give the experiment design a name (i.e. AC) and then use the **Browse** function to locate the Tutorial Groups.csv file.



The 'Group runs by' drop down will update to reflect the possible fields in the csv file that you can use to 'group' your runs by: Conditions, Date of collection or Location.



Select Conditions and then click **Next**.

**Note:** you can create additional experiment designs following the completion of automatic processing.


The Protein Quantitation dialog opens displaying the default method, **Relative Quantitation using Hi-N** which uses up to 3 peptides per protein to compute the relative amount of each 'identified' protein.




Use the drop down to reveal the alternative methods for protein quantitation

Select the Default option **Relative Quantitation using Hi-N**, with **Use protein grouping** ticked then click **Finish**.

The Alignment process starts with the automatic selection of A2 as the alignment reference



Once Alignment completes Peak Picking commences




When Processing completes, depending on what stages you selected to perform, the Wizard displays (on the right) what stage the workflow will open at.

In this example, as the data is DDA it is going to open at **Identify Peptides** if all possible processing steps are selected.

You can either:

- Continue with the analysis, as the Processing dialog is not displaying any warnings, and perform **Identify Peptides**. In which case you can go to page 41
- Open the analysis at Identify Peptides and immediately return to the Review Alignment stage by clicking on it in the Workflow to review alignment quality (page 14)
- Alternatively, you can **Close** the dialog. This will **not** move you to a later stage in the workflow but instead will allow you to review the 'Run Chromatography' then navigate systematically through all the stages yourself or jump from stage to stage.




**Note:** if processing fails to complete successfully there are a number of suggested strategies you can use to proceed with your analysis. These are outlined in Appendix 2 (page 76)

The next stage in this document gives you a graphical view of how to proceed with your analysis following automatic processing.

## Stage 2B: After Automatic Processing

When Processing completes, depending on what stages you selected to perform, the Automatic Processing Dialog displays what stage the workflow will open at.




The flow chart gives you an overview of the various steps you can take to proceed with your automatically processed data.

In this example, as the data is DDA it is going to open at Identify Peptides so you can either:

- Continue with the analysis, as the Processing dialog is not displaying any warnings, and perform **Identify Peptides**. In which case you can go to page 41
- Open the analysis at Identify Peptides and return to the Review Alignment stage by clicking on it in the Workflow (page 14).

## Stage 3: Licensing

This stage in the analysis workflow will **only** appear if you are using 'Unlicensed' data files to evaluate the software and have no dongle attached.



If you have performed an analysis using Automatic Processing without a valid dongle or do not have the appropriate code to licence your runs, and you close Progenesis QI for proteomics you will be warned that the analysis will be lost.

**QIP Progenesis QI.p Tutorial for DDA - Progenesis QI for proteomics**

File

Import Data Licensing Review Alignment Filtering Experiment Design Setup Review Peak Picking Peptide Ion Statistics Identify Peptides QC Metrics Refine Identifications Review Proteins Protein Statistics Report

nonlinear  
A Waters Company

**Dongle License Runs**

This installation is currently restricted to analyse licensed runs only.

To license your runs, you need an evaluation or lease licence code which can be obtained from a sales representative.

Once licensed, your runs can be analysed on any installation of the software. The licence is automatically included when archiving an experiment.

If your runs have been licensed on another computer, [click here](#) to make the licences available on this computer.

If you have one, you can [open a licence file](#) to install.

If you have just installed a dongle, [click here](#).

Run name	Licence state	License this run
C:\Users\andy.borthwick\Documents\Customer Data\Progenesis QI.p v3.0 Tutorials and Demo Suites\Proge...	Unlicensed	<input checked="" type="checkbox"/>
C:\Users\andy.borthwick\Documents\Customer Data\Progenesis QI.p v3.0 Tutorials and Demo Suites\Proge...	Unlicensed	<input checked="" type="checkbox"/>
C:\Users\andy.borthwick\Documents\Customer Data\Progenesis QI.p v3.0 Tutorials and Demo Suites\Proge...	Unlicensed	<input checked="" type="checkbox"/>
C:\Users\andy.borthwick\Documents\Customer Data\Progenesis QI.p v3.0 Tutorials and Demo Suites\Proge...	Unlicensed	<input checked="" type="checkbox"/>
C:\Users\andy.borthwick\Documents\Customer Data\Progenesis QI.p v3.0 Tutorials and Demo Suites\Proge...	Unlicensed	<input checked="" type="checkbox"/>
C:\Users\andy.borthwick\Documents\Customer Data\Progenesis QI.p v3.0 Tutorials and Demo Suites\Proge...	Unlicensed	<input checked="" type="checkbox"/>

Run licence code:  **Use Licence Code**

**Section Complete** →

**QIP Unable to save experiment**

You cannot save the analysis without a valid license. If you close now your analysis will be lost.

OK Cancel

**Note:** although the analysis will be lost (Alignment Peak Picking etc) the experiment file and all the 'loaded pre-processed runs are retained'. This allows you to reopen the experiment, when you have a licenced code and/or dongle available, and redo the analysis steps without having to reload the data.

For more details on managing the licencing of your data refer to Appendix 4 (page 79)

**If you are using the tutorial archive, this page will not appear as the data files are licensed.**

## Stage 4: Review Alignment

At this stage Progenesis QI for proteomics Alignment opens displaying the alignment of the runs to the Reference run (A2).

### Layout of Alignment

To familiarize you with Progenesis QI for proteomics Alignment, this section describes the various graphical views used in the alignment of the LC-MS runs.

To setup the display so that it looks similar to the one above:

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/> <input type="checkbox"/>	416	98.9%
A2	<input checked="" type="checkbox"/> <input type="checkbox"/>	Ref	
A3	<input checked="" type="checkbox"/> <input type="checkbox"/>	413	99.2%
C1	<input checked="" type="checkbox"/> <input type="checkbox"/>	223	79.9%
C2	<input checked="" type="checkbox"/> <input type="checkbox"/>	227	75.8%
C3	<input checked="" type="checkbox"/> <input type="checkbox"/>	243	79.1%

- In the Run table click on Run C1 to make it current. You will now be looking at the alignment of C1 to A2 in the unaligned view. Now drag out an area to review on the **Ion intensity map**. The other 3 views will update to reflect the new focus.

**Reference Run (Magenta)**

**Current Run (Green)**

**Table of Alignment Vectors and Scores**

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/> <input type="checkbox"/>	416	98.9%
A2	<input checked="" type="checkbox"/> <input type="checkbox"/>	Ref	
A3	<input checked="" type="checkbox"/> <input type="checkbox"/>	413	99.2%
C1	<input checked="" type="checkbox"/> <input type="checkbox"/>	223	79.9%
C2	<input checked="" type="checkbox"/> <input type="checkbox"/>	227	75.8%
C3	<input checked="" type="checkbox"/> <input type="checkbox"/>	243	79.1%

**Vector Editing (Window A):** is the main alignment area and displays the area defined by the current focus rectangle shown in Window C. The current run is displayed in green and the reference run is displayed in magenta. Here is where you can review in detail the vectors and also place the manual alignment vectors when required.

**Transition (Window B):** uses an **alpha blend** to animate between the current and reference runs. Before the runs are aligned, the peptide ions appear to move up and down. Once correctly aligned, they will appear to pulse. During the process of adding vectors, this view can be used to zoom in and also navigate thus helping with accurate placement of manual vectors.

**Ion intensity map (Window C):** displays the alignment quality across retention time and m/z. A blue box indicates the current focus area.

**Total ion chromatogram (Window D):** shows the intensity of the runs over time, with peaks corresponding to the retention times of the ions.

Alignment quality: ■ Good ■ OK ■ Needs review

Section Complete

**Vector Editing (Window A):** is the main alignment area and displays the area defined by the current focus rectangle shown in Window C. The current run is displayed in green and the reference run is displayed in magenta. Here is where you can review in detail the vectors and also place the manual alignment vectors when required.


**Transition (Window B):** uses an **alpha blend** to animate between the current and reference runs. Before the runs are aligned, the peptide ions appear to move up and down. Once correctly aligned, they will appear to pulse. During the process of adding vectors, this view can be used to zoom in and also navigate thus helping with accurate placement of manual vectors.

**Ion Intensity Map (Window C):** shows the **focus** for the other windows. When you click on the view the blue rectangle will move to the selected area. The focus can be moved systematically across the view using the cursor keys. The focus area size can be altered by clicking and dragging out a new area with the mouse. This view also provides a visual quality metric for the Alignment of the runs, from green through yellow to red.

**Total Ion Chromatogram (Window D):** shows the current **total ion** chromatogram (green) overlaid on the Reference chromatogram (magenta). As the peptide ions are aligned in the **Vector Editing** view the chromatograms become aligned. The retention time range displayed is the vertical dimension of the Focus Grid currently displayed in the **Ion Intensity map** (Window C).

## Reviewing quality of alignment vectors

After **Automatic alignment** the number of vectors and Quality Scores will be updated on the **Runs** panel and the vectors will appear (in blue) on the **Vector editing** view.



If the alignment has worked well, the Ion Intensity Map will appear predominantly green and a score of greater than 80% will be reported in the table. Also, when **Show Aligned** is selected, in Window A (vector editing) vector length should appear minimal and in Window B (Transition) will show peptide ions pulsing slightly but not moving up and down.



**Note:** you can use the icon to the right of the panel name to expand or contract each view.

## Reviewing Quality of Alignment



At this point the quality metric, overlaid on the Ion Intensity Map as coloured squares, acts as a guide drawing your attention to areas of the alignment. These range from Good (Green) through OK (Yellow) to Needs review (Red). When reviewing individual squares double click on a coloured square to set the focus.

Three example squares are examined here.



For a 'green' square, the majority of the data appears overlapped (black) indicating good alignment. When viewed in the Transition view the data appears to pulse.



For a 'yellow' square, some of the data appears overlapped (black) indicating OK alignment. When viewed in the Transition view some of the data appears to pulse.



For a 'red' square little of the data appears overlapped (black) indicating questionable alignment. When viewed in the Transition view little data appears to pulse.




**Note:** the coloured metric **should be used as a guide**. In cases where there are a few 'isolated' red squares this this can also be indicative of 'real' differences between the two runs being aligned and should be considered when examining the overall score and surrounding squares in the current alignment.

The weighted average of the individual squares gives the overall percentage score for each run alignment.



**Note:** a marked area of red combined with a low score clearly indicates a 'misalignment' and may require some manual intervention (see Appendix 5, page 80).



Having performed the analysis automatically, in the course of reviewing the quality of alignment you decide that the alignment requires editing then you can unlock the analysis.

To do this click on the **This run is protected from editing** link above the 'greyed out' **Align runs automatically** button.

**Review Alignment**  
Sample ions are aligned to compensate for drifts in retention time between runs.

[Learn about the visualisations shown here](#)

**1 Align retention times automatically**  
For maximum reproducibility, the software can automatically align your runs.

This run is protected from editing

**2 Review the alignment**  
Using the quality control measures, review and edit the runs' alignment:

- Order the runs by alignment score and start by selecting the first run
- Within each run, inspect and edit any areas rated as Needs Review

[Learn about the review and editing process](#)

**QI Protected from editing**

**Delete existing analysis?**

If you change the alignment, it will invalidate the current analysis including peptide ion pattern and IDs, editing, and tags.

If you want to keep these, you should archive this experiment before changing the alignment and moving forward.


**Note:** as you click on the link, you will be warned that you are discarding the current analysis (all steps beyond alignment).

Details on editing alignment are described in Appendix 5 (page 80)

**The alignment quality of this tutorial data set does not require any manual intervention.**

If you have chosen to discard the current automatic analysis and have been exploring the alignment of one or more of the runs using manual vectors (for the purposes of this tutorial) make sure you first remove all manual vectors and then re-perform the Automatic alignment.

To do this for C1, first select Remove 'All vectors in the whole run' and then click Align runs automatically.



Having re-performed the Automatic alignment the Vectors and scores will appear as shown above.



Since you have unlocked your analysis you must now perform the Peak Picking that was originally performed as part of the **Automatic processing**.

To move to Peak Picking click on **Filtering** on the Workflow or **Section Complete** (bottom right).

## Stage 5A: Filtering

The Peak picking Parameters dialog opens if Picking has not been performed. If it has been performed move to section 5B.

### Peak Picking Parameters



The Peak Picking Parameters dialog opens, showing all the runs in the current experiment and a tick against each run. This is the default setting, where the peak picking algorithm uses information from all of the runs to contribute to the pattern of peptide ion outlines.

**Tip:** It may be appropriate **only** to pick peaks that are present in a limited number of your runs. In which case un-tick the runs that you do **NOT** want to contribute to the peptide ion detection pattern. This may be important when one or more of the 'replicate' runs appear noisy due to non-optimal chromatography or sample handling.

**Note:** peptide ions outlines will be added to 'un-ticked' runs; although these runs will **not** contribute to the peak picking pattern.

**Tip:** depending on run quality, a suggested minimum number of ticked runs should include at least one replicate of each experimental condition.

The sensitivity of the peak picking can be controlled by adjusting settings under the **Peak picking limits** tab.

**Note:** for all 3 **Sensitivity** modes a Chromatographic peak width (Retention time window) for the peaks **can** be set by applying a minimum retention window or peak width in minutes.

**Chromatographic peak width**

The chromatographic peak width gives the length of time over which an ion has eluted. If you set a minimum peak width, any ion that has eluted over a shorter period will be rejected.

Apply a minimum peak width

Minimum width:  minutes

**For the runs in this user guide we will use the default settings for the Automatic method and NOT apply a minimum peak width.**

The third tab allows you to set the **maximum charge** of the peptide ions, which will be detected. The default setting is a charge state of 20. For this example, leave this set as default.

**Note:** you can either, leave this set as default and remove the high charge state peptide ions at the Filtering stage or you can choose not to detect them in the first place by reducing the charge state threshold here.

Finally, you can set **Retention time limits** for the detection. Default values are displayed in the unticked boxes and detection will use these values whenever the box(s) are unticked. Enter values of 10 and 75 min and tick the boxes as shown below.

QI/P Peak Picking Parameters

Runs for peak picking | Peak picking limits | **Maximum charge** | Retention time limits

**Maximum allowable charge**

You can set the maximum charge of ions to be detected. Ions with a charge greater than this value will be rejected.

Maximum ion charge:

Start peak picking | Cancel

QI/P Peak Picking Parameters

Runs for peak picking | Peak picking limits | Maximum charge | **Retention time limits**

**Retention time limits**

You can set the minimum and maximum retention time for peak picking. Ions that elute before or after these values will be ignored.

Ignore ions before:  minutes

Ignore ions after:  minutes

Start peak picking | Cancel

Click **Start Peak Picking** to start the detection process.

On completion of detection, the Filtering stage will open displaying the number of peptide ions detected, in this example there are 14473 peptide ions.

If required you can remove peptide ions based on position, charge state, number of isotopes or combinations of these peptide ion properties.

As an example, we will filter the peptide ions based on 'charge state'.

When **With charge** is selected the number of peptide ions present at each charge state is displayed, these can be selected accordingly. In this case we will retain peptide ions with a charge state of 2 to 7.

Area limits, charge state and number of isotopes can be combined to refine the peptide ion selection.

**Tip:** when filtering on only one property of the peptide ion i.e. charge state, make sure you have 'collapsed' the other filters (see right)

Hence all peptide ions with a charge state of 1 or 8 and above will appear **red** on the main view as you hold the cursor over the appropriate delete button.

The screenshot displays the Progenesis QI software interface during the filtering stage. The left sidebar shows a list of charge states with checkboxes: Charge 5 (465 peptide ions), Charge 6 (79 peptide ions), Charge 7 (37 peptide ions), Charge 8 (21 peptide ions), Charge 9 (16 peptide ions), Charge 10 (14 peptide ions), Charge 11 (8 peptide ions), Charge 12 (9 peptide ions), Charge 13 (5 peptide ions), Charge 14 (4 peptide ions), Charge 15 (3 peptide ions), Charge 16 (5 peptide ions), Charge 17 (7 peptide ions), Charge 18 (5 peptide ions), Charge 19 (2 peptide ions), and Charge 20 (1 peptide ion). Below this list, there is a 'Number of isotopes' dropdown, a 'Delete 13559 matching peptide ions' button, and a highlighted 'Delete 914 non-matching peptide ions' button. A note indicates '(14473 peptide ions in total)'. Below the note, there is a 'Normalisation' section with a note: 'Note that all remaining peptide ions will be used in the normalisation calculation. If deleting a substantial number of peptide ions, you should review the normalisation afterwards.' and a 'Review normalisation >>' button. The main plot shows a dense cloud of points representing peptide ions, with retention time (min) on the y-axis (ranging from 20.8 to 99.2) and m/z on the x-axis (ranging from 612.5 to 1887.5). The 'Undo' button is highlighted with a red box. A 'Section Complete' button is visible in the bottom right corner.

To remove these peptide ions press **Delete 914 Non Matching Peptide ions**.

**Note:** if you have deleted peptide ions unintentionally you can use the **Undo** button to bring them back, however, when you move to the next section you will lose the capacity to undo the filter.


In this example do not undo the charge filter you performed as described above.

Before moving on from filtering you can review the normalisation of the data.


**Tip:** When you have reached the filtering stage, it is good practice to close the experiment and save an archive. This can be used to restore the unfiltered state if the filtering you have performed is too 'stringent'.

## Stage 5B: Reviewing Normalisation

Review normalisation is accessed from the button at the bottom left corner of the filtering page.



If you have filtered out a number of peptide ions from the original detection pattern then the normalisation will update.




The **Review Normalisation** page will open displaying plots for the normalisation of all the peptide ions on each run.

This page in the workflow **does not** allow you to alter the Normalisation of your data but provides you with individual views for each run showing the data points used in the calculation of the normalisation factor for the run.

Alternatively, if you do not believe normalisation is necessary, you can opt to 'Don't use any normalisation' for the rest of the analysis (Normalisation Method tab).


Normalisation factors are reported in the table to the left of the plots.



### Calculation of Normalisation Factor:


Progenesis QI for proteomics will automatically select one of the runs that is 'least different' from all the other runs in the data set to be the 'Normalising reference'. The run used, is shown above the table of Normalisation factors.

For each sample run, each blue dot shows the log of the abundance ratio for a different peptide ion (normalisation target abundance/run abundance).




The details for individual peptide ions can be viewed as you hold the cursor over the dots on the plot.

On the graph the peptide ions are shown ordered by ascending mean abundance. The normalisation factor is then calculated by finding the mean of the log abundance ratios of the peptide ions that fall within the 'robust estimated limits' (dotted red lines). Peptide ions outside these limits are considered to be outliers and therefore will not affect the normalisation.




**Note:** if you do **not** wish to work with normalised data then Select '**Don't use any normalisation**'.



**Note:** once you have identified a set of proteins, you can then apply the **Normalise to a set of housekeeping proteins** by using this option to locate and select the peptide ions.

For this experiment, you should leave the **Normalise to all peptide ions** option selected.

Now return to filtering by clicking on the button on the bottom left of the screen



For this example, we **DO NOT** do any additional Filtering beyond the Charge state filter you applied on page 21, so click on **Section complete**.

**Note:** if you do any extra filtering then **Normalisation recalculates** as you move to the next stage in the Workflow.




## Stage 6: Experiment Design Setup for Analysed Runs

At this stage in the workflow you can setup one or more experimental designs for your data.

There are two basic types of experimental designs:

**Between-subject design:** here samples from any given subject appear in only one condition. (i.e. control versus various drug treatments). The ANOVA calculation assumes that the conditions are independent and applies the statistical test that assumes the means of the conditions are equal.




**Within-subject design:** here samples have been taken from a given subject under different conditions (i.e. the same subject has been sampled over a period of time or after one or more treatments). Here a standard ANOVA is not appropriate as the data violates the ANOVA assumption of independence. Therefore by using a *repeated measures* ANOVA, individual differences can be eliminated or reduced as a source of between condition differences. This within-subject design can be thought of as an extension of the paired samples t-test, including comparison between more than two repeated measures.


Additional information on how to apply the Within-subject Design is in Appendix 6 page 84

This experiment contains 2 conditions: A and C and uses the **Between-subject design** to group the analysed runs to reflect the Biological conditions in the original study.

To create a new **Between-subject Design** move the cursor over this option and click to open the dialog.




Give the new experimental design a name and then click **Create design**.



## To create a new condition


1. Select the runs for the condition by clicking on the required icon in the **Runs** panel, as shown.
2. Press the 'black triangle' next to the **Add Selected Runs to Condition** button on the main toolbar.
3. Select **Add to new condition...** from the drop down menu.
4. A new condition will appear in the **Conditions** panel on the left.
5. Rename the condition (e.g. C) by over typing the default name.
6. Repeat steps 1 to 5 until all the required runs are grouped into conditions.

An alternative way to handling the grouping of this set and other larger (and more complex) experimental designs is to make use of **sample tracking information** that has been stored in a spreadsheet at the time of sample collection and/or preparation.



For this example there is a **Tutorial Groups.csv** file available in the Experiment Archive you restored at the beginning of this tutorial exercise.


To use this approach select the **Import design from file** option from the New Experiment Design dialog. Then locate the Tutorial Groups file and select what to **Group by**, for example: **Conditions**.



When **Create design** is pressed the new tab refreshes to allow you to adjust the conditions.

**Note:** currently this second method is the only method of creating an Experiment Design that can be applied when using the Automatic Analysis process. Additional Designs can be applied post Automatic processing.

Use **Delete** on the Conditions panel to remove conditions that are not required in this particular design.



**Note:** On deleting a condition the runs will reappear in the Runs window.

**Note:** both designs are available as separate tabs.

To move to the next stage in the workflow, **Review Peak Picking**, click **Section Complete**.

## Stage 7: Review Peak Picking and editing of results

The purpose of this stage in the Workflow is to review the list of peptide ions using the visual tools provided and edit peptide ions if required.

The review stage has 4 display modes: 1D, 2D, 3D and Peptide ion Details controlled by the tabs on the bottom left of the display and the expander bar to the right of the table. Each display has multiple views to allow comparative exploration of the detected peptide ions on the aligned LC-MS runs.

### Exploring analysed data using the Data displays

**Window A:** shows the list of peptide ions ranked by the p value for the one way **Anova** using the current grouping.

**Note:** a value of 'Infinity' in the **Fold** column indicates 'Presence/Absence'

The screenshot displays the Progenesis QI software interface for the 'Review Peak Picking' stage. The main window is titled 'Progenesis QI.p Tutorial for DDA - Progenesis QI for proteomics'. The interface is divided into several sections:


- Menu Bar:** Includes 'File', 'Review Alignment', 'Filtering', 'Experiment Design Setup', 'Review Peak Picking', 'Peptide Ion Statistics', 'Identify Peptides', 'QC Metrics', 'Refine Identifications', 'Review Proteins', 'Protein Statistics', and 'Report'.
- Review Peak Picking Panel:**
  - Buttons: 'No filter applied', 'Create...'
  - Table of peptide ions:
- Mass spectrum (26.424 min):** A plot showing intensity versus m/z, with a prominent peak at m/z 1099.5687. Labeled 'B'.
- Chromatogram (m/z = 1099.5687):** A plot showing intensity versus retention time (min), with a peak at 26.424 min. Labeled 'C'.
- 1D Display:** A section with buttons for '1D Display', '2D Display', and '3D Display'. Labeled 'D'.
- 2D Display:** A scatter plot showing intensity versus retention time (min) and m/z. Labeled 'E'.
- Run Section:** Includes 'Expression profile', 'Run:' dropdown, 'Aggregate' dropdown, and buttons for 'Select', 'Edit', 'Split', 'Merge', 'Add', 'Delete', 'Undo', and 'Redo'. A 'Charge state' legend is also present.
- Bottom Bar:** Includes 'Experiment design: AC' and a 'Section Complete' indicator.

#	Anova (p)	q Value	Fold	Tag	Not
5770	1.14E-10	1.46E-08	Infinity		
14169	1.16E-10	1.46E-08	Infinity		
2725	1.23E-10	1.46E-08	Infinity		
1656	1.25E-10	1.46E-08	Infinity		
1613	1.43E-10	1.59E-08	Infinity		
14174	2.06E-10	1.99E-08	Infinity		
2635	2.06E-10	1.99E-08	1.74E+07		
11235	2.08E-10	1.99E-08	Infinity		
12123	2.14E-10	1.99E-08	Infinity		
641	2.36E-10	2.1E-08	Infinity		


**Note:** by default all the peptide ions are included in the selection for the next section of the analysis.

## The 1D Display

**Window B:** displays the Mass spectrum for the current peptide ion on the selected Run (in window D). Hold the cursor over the peak to indicate m/z and intensity



**Window C:** displays the Chromatogram for the current peptide ion on the selected Run (in window D). Hold the cursor over the peak to indicate Retention time and intensity





**Window D:** displays the details of the currently selected run. By default the selected run is an Aggregate of all the aligned runs.

Details of individual runs can be viewed by using the 'Run' link and selecting the run you wish to view.

The peptide ion editing tools are located in this window (see page 32 for functional explanation).


Clicking on the Expression Profile tab in Window D shows the comparative behaviour of the peptide ion across the various biological groups based on group average normalised volume. The error bars show +/- 3 standard errors.




**Window E:** shows where the current peptide ion is located on the LC-MS run by means of the cross hairs.

To change the current location, click on the image of the run (**note:** the retention time and m/z values update as you move the cursor around this view).

**Note:** doing this updates the focus of all the other windows.




You can also drag out a larger area on this view that will refocus the other windows.



### The 2D Display


Windows A, D and E: perform the same functions across all 4 display modes.

In the 2D Display mode, Window B displays a montage of the current peptide ion across all the aligned LC-MS runs.



The appearance of the Montage (window B) is controlled by the panel on the bottom left of the display.


Using the the various views in the 2D display one can examine the peptide ion detection in detail to validate the correct detection of even fully overlapping peptide ions as shown below.



**Note:** Use **Go To Location** tool, in the Run window, to find peptide ion of interest then click on the monoisotopic envelope to refocus view.

### The 3D Display

Window B changes into a 3D view by selecting the 3D Montage tab on the bottom left of the display.




The views can be set to **Rotate** automatically or you can rotate them manually by clicking and dragging them with the mouse.


## Editing of peptide ions in the View Results stage

As an example of using the editing tools which are located on the left of the **Run** view, we will remove and add back the 'monoisotopic peak' for the detected peptide ion selected below. A peptide ion can be selected from the 'Peptide ions' list or located using the various views.


1. Locate the peptide ion at approx 1005.4 m/z and 31.7 min using the **Go To Location** tool (at the intersection of the m/z and RT axis).



2. Select the **Edit** tool and click on the peptide ion to reveal the 'edit handles'




3. Click on the 'minus' handle over the monoisotopic peak to remove it.




4. Click outside the boundary of the peptide ion to update the view.



- To add a peak to an existing peptide ion, ensure that **Edit** is selected then click inside the peptide ion to reveal the handles.



- Click on the 'plus' handle on the peak to add it.
- Then click outside the peptide ion to update the view.
- Note:** If you are not satisfied with the editing use the **Undo** button and retry.
- Note:** that a tag is automatically added to the edited peptide ion in the table and the peptide ions id number is changed to the next available one at the end of the list.



**Note:** use **Undo** to reverse this editing process, this restores the original # and removes the **Edit** tag.


Also to remove the Edit tag right click on the table, select Edit tags, and delete the Edit Tag from the list.

The other tools: split, merge, add and delete behave in a similar fashion and their use can be combined to achieve the desired results.

### Selecting and tagging peptide ions for Peptide Ion Statistics

There are a number of ways to 'refine' your 'Ranked List' of analysed peptide ions before examining them with the Statistical tools in **Peptide Ion Statistics**. These make use of simple 'Selection' and 'Tagging' tools that can be applied to the various groupings created in Stage 6 (page 25). An example is described below.

First expand the 'Peptide ions' table to show all the details by double clicking on the 'Splitter Control' to the right of the Review Peptide ions table.




Then order on **Abundance** and select all peptide ions with an Abundance of 1E+05 and greater, (the exact number is not important).

#	Anova (p)	q Value	Fold	Tag	Notes	Highest Mean	Lowest Mean	m/z	z	Mass	RT (mins)	RT window (mins)	Abundance	Intensity	Max CV (%)	MS/MS	Protein	Peptide Score	Peptide Mod
1	0.414	0.0823	1.05		Add a note...	A	C	805.4410	3	2413.301	55.594	7.27	1.24E+08	1.08E+08	7.22	152			
17	0.183	0.0412	1.1			C	A	1207.6552	2	2413.296	55.594	4.64	9.37E+07	5.04E+07	7.62	102			
39	5.49E-08	4.05E-07	1.16E+05			C	A	1100.5857	3	3298.735	45.749	3.15	8.53E+07	7.87E+07	15.9	19			
61	2.34E-06	5.8E-06	1.6E+03			A	C	1176.2272	3	3525.660	48.790	8.8	6.83E+07	2.22E+07	23.5	55			
9	3.43E-08	3.05E-07	1.35E+03			A	C	656.8613	2	1311.708	44.596	4.21	6.1E+07	1.17E+08	8.98	39			
19	1.87E-07	9.36E-07	805			C	A	988.9848	2	1975.955	51.484	4.82	5.3E+07	9.19E+07	14.5	48			
10	6.3E-08	4.29E-07	134			C	A	663.8693	2	1325.724	47.557	4.77	5.18E+07	1.69E+08	6.45	67			
23	3.4E-07	1.41E-06	6.37E+03			C	A	900.9713	2	1799.928	40.004	2.87	4.71E+07	8.16E+07	24.6	15			
54	3.63E-06	8.06E-06	3.25E+03			A	C	1061.0071	2	2120.000	54.334	6.1	4.2E+07	2.52E+07	26.6	64			
38	5.04E-06	1.05E-05	1.28E+03			A	C	997.4478	2	1992.881	32.637	2.53	4.19E+07	3.81E+07	34.1	29			

With the peptide ions still highlighted right click on them and select 'New Tag'.

#	Anova (p)	q Value	Fold	Tag	Notes	Highest Mean	Lowest Mean	m/z	z	Mass	RT (mins)	RT window (mins)	Abundance	Intensity	Max CV (%)	MS/MS	Protein	Peptide Score	P
2286	0.535	0.101	1.15		Add a note...	A	C	705.3259	2	1408.637	47.144	2.83	1E+05	3.41E+05	44.2	4			
3496	0.162	0.0363	1.64		Add a note...	C	A	925.7347	4	3698.910	47.282	0.86	1E+05	2.48E+05	44.2	1			
6490	0.024	0.0082	11.7		Add a note...	C	A	957.7371	4	3826.919	54.884	1.05	1E+05	1.51E+05	100	1			
5200	0.0414	0.0124	148		Add a note...	A	C	1061.0065	2	2119.999	49.045	0.644	1E+05	1.44E+05	126	5			
1330	0.00682	0.00316	2.05		Add a note...	C	A	42.014	0.935		42.014	0.935	1E+05	8.97E+05	21.5	10			
1439	0.0267	0.00888	2.1			A	C	39.518	1.05		39.518	1.05	9.99E+04	4.53E+05	26.2	9			
4175	0.0169	0.00634	392			C	A	35.643	0.595		35.643	0.595	9.98E+04	2.69E+05	173	2			
1244	0.00648	0.00304	2.42			A	C	42.725	0.817		42.725	0.817	9.98E+04	5.23E+05	22.4	7			
4973	0.0536	0.015	126			C	A	54.032	0.831		54.032	0.831	9.98E+04	1.96E+05	87.7	0			
2085	0.0684	0.0181	3.29			C	A	943.1200	3	2826.338	56.463	1.35	9.97E+04	4.53E+05	60.8	0			

Give the Tag a name. i.e. **'Most abundant'**.




On clicking **OK** the Tag is added to the peptide ions highlighted in the table (signified by a coloured square).

Now right click on any peptide ion in the table and select **Quick Tags** this will offer you a number of standard tag options. Select **Anova p-value....** Then set the threshold as required and adjust default name as required and click **Create Tag**.


#	Anova (p)	q Value	Fold	Tag	Nc
3496	0.162	0.0363	1.64	■	
6490	0.024	0.0082	11.7	■	
5200	0.0414	0.0124	148	■	
1330	0.00682	0.00316	2.05	■	
1439	0.0267	0.00888	2.1		
4175	0.0169	0.00634	392		
1244	0.00648	0.00304	2.42		
4973	0.0536	0.015	126		
2085	0.0684	0.0181	3.29		
1381	3.15E-06	7.27E-06	7.84		

#	Anova (p)	q Value	Fold	Tag	Notes	Highest Mean
3496	0.162	0.0363	1.64	■		C
6490	0.024	0.0082	11.7	■		C
5200	0.0414	0.0124	148	■		A
1330	0.00682	0.00316	2.05	■		C
1439	0.0267	0.00888	2.1			A
4175	0.0169	0.00634	392			
1244	0.00648	0.00304	2.42			
4973	0.0536	0.015	126			
2085	0.0684	0.0181	3.29			
1381	3.15E-06	7.27E-06	7.84			




Once this tag appears against peptide ions in the table right click on the table again and create another Quick Tag, this time for peptide ions with a **Max fold change ≥ 2**

The table now displays peptide ions with multiple tags.






The tags can be used to quickly focus the table on those peptide ions that display similar properties.

For example: to focus the table on displaying those peptide ions that have an **Anova p-value ≤ 0.05** click on **Create** on the filter panel above the table.



#	Anova (p)	q Value	Fold	Tag	Not
1613	1.43E-10	1.59E-08	Infinity	■	
927	0.284	0.0583	1.25	■	
489	7.52E-05	9.24E-05	11.5	■	
2050	0.0491	0.014	8.71	■	
1067	0.000316	0.000285	7.79	■	
1743	0.0144	0.00562	3.29	■	
403	0.865	0.152	1.02	■	
2529	3.94E-05	5.52E-05	155	■	
426	0.0489	0.014	1.86	■	
469	0.0331	0.0105	1.48	■	


Drag the tag on to the panel **Show peptide ions that have all of these tags** and press OK.




Now order the current peptide ions in the table by the **Highest mean** so that all the peptide ions showing the highest mean for **Condition C** are at the top of the list.

Then highlight all the peptide ions with the highest mean for **condition C** and create a new Tag for them.

#	Anova (p)	q Value	Fold	Tag	Notes	Highest Mean
387	2.63E-05	3.99E-05	21.4		Add a note...	C
12372	0.0337	0.0106	169		Add a note...	C
2946	0.00336	0.00181	5.4		Add a note...	C
2953	0.00923	0.004	1.02E+03		Add a note...	C
12376	0.0107	0.00449	668		Add a note...	C
9122	0.0127	0.00512	354			A
14212	0.000258	0.000244	Infinity			A
9120	0.0251	0.00849	86.6			A
14407	0.0266	0.00886	2.77			A
9307	0.0062	0.00293	6.75			A




#	Anova (p)	q Value	Fold	Tag	Notes	Highest Mean
387	2.63E-05	3.99E-05	21.4		Add a note...	C
12372	0.0337	0.0106	169		Add a note...	C
2946	0.00336	0.00181	5.4		Add a note...	C
2953	0.00923	0.004	1.02E+03		Add a note...	C
12376	0.0107	0.00449	668		Add a note...	C
9122	0.0127	0.00512	354			A
14212	0.000258	0.000244	Infinity			A
9120	0.0251	0.00849	86.6			A
14407	0.0266	0.00886	2.77			A
9307	0.0062	0.00293	6.75			A




Create a tag for them called **Significantly up in C**, tagging 3946 peptide ions.

Finally view the tags you have just created by clicking on Edit in the Tag filter panel, above the table.



Make sure that only the tag for the **Most abundant** peptide ions is shown and press **OK**.

To move to the next stage in the workflow, Peptide Statistics, click **Section Complete**.



## Stage 8: Peptide Ion Statistics on selected peptide ions


The user guide now describes the functionality of the Multivariate Statistics.

Peptide Statistics opens calculating the Principal Components Analysis (PCA) for the active 'tag' in this case the **Most abundant** peptide ions.



As an example we will start by examining the behaviour of the **Most abundant** peptide ions from the previous stage, **Review Peak Picking**.

The statistical analysis of the selected data is presented to you in the form of interactive graphical representation of answers to questions asked of the analysed data.



**Note:** the LC-MS runs (samples) are displayed as solid coloured circles on the plot. To identify the runs, a tooltip is displayed when the cursor is held over each circle.


### Principal Component Analysis (PCA)

In **Peptide Stats** the first statistically based question asked of the data takes the form of a Quality Control assessment:

*Are there any outliers in my data? And does my data cluster according to my experimental conditions?*

It answers this question by:

*'Using Principal Components Analysis (PCA) to produce a simplified graphical representation of your multidimensional data'.*



PCA can be used to determine whether there are any outliers in the data and also look at how well the samples group. The groupings observed on the 2D PCA plot can be compared to your experimental conditions and conclusions made regarding possible outliers in your data. Selecting peptide ions in the table will highlight the peptide ions on the 'Biplot' and their expression profiles will appear in the lower panel.

**Note:** the Table in the Stats view contains additional columns:

**q value:** tells us the expected proportion of false positives if that peptide ion's p-value is chosen as the significance threshold.

**Power:** can be defined as the probability of finding a real difference if it exists. 80% or 0.8 is considered an acceptable value for power. The Power Analysis is performed independently for each peptide ion, using the expression variance, sample size and difference between the means.

Also, for a given power of 80% we can determine how many samples are required to ensure we find a difference if it actually exists.


**Note:** Power analysis is discussed in Appendix 7 (page 86)

## Correlation Analysis

Use the tags, created in Review Peak Picking, to filter the peptide ions in the table.


We are going to explore the Correlation Analysis for all the peptide ions that were tagged at the view results stage with an **Anova p-value  $\leq 0.05$** .

On pressing OK the PCA will recalculate using these peptide ions, you can (to save time) stop this calculation by pressing **Cancel calculation** and then set up Correlation Analysis for the 8174 peptide ions.




To set up the **Correlation Analysis** using this filtered data set click on **Ask another question** (above the table)

A selection of 3 tools will appear in the form of questions.



Select the second option to explore 'peptide ion correlation based on similarity of expression profiles'




This time the statistically based question(s) being asked is:

*'Group my (selected) peptide ions according to how similar their expression profiles are'*

The question is answered by:


*'Using Correlation analysis to evaluate the relationships between the (selected) peptide ions' expression profiles'*

The answer is displayed graphically in the form of an interactive dendrogram where the vertical distance, between each peptide ion, can be taken as indicative of how similar the expression profiles of each cluster of peptide ions are to each other.




**Correlation Analysis** enables the grouping of peptide ions together according to how similar their expression profiles are.

To highlight all the peptide ions demonstrating **Increased expression in the C** group click on a 'node' for a branch of the Dendrogram (as shown above). As before, right click on the highlighted peptide ions in the table and create a Tag for these peptide ions (Up regulated in C).




Also create a tag for those peptide ions showing **Increased expression in A** by first clicking on the other 'main' node then right click on the highlighted peptide ions in the table and creating the New tag (Up regulated in A).



**Comment:** When you review the tags using **Edit** you can see that the Magenta and Yellow tags have been assigned to the same number of peptide ions. This shows how tabulated information about peptide ions can be used alongside interactive graphical plots of multivariate statistical analysis to explore your data.

**Note:** two groups is a special case, for more groups this will not be the case, and additional tagging will be required.





**Tip:** when reviewing the tags (see above) if you are not applying a new filter then use the **Cancel** button to return to the main view, this prevents unnecessary recalculation of your data.

To move to the next stage in the workflow, **Identify Peptides**, click **Section Complete**.

## Stage 9: Identify peptides

Progenesis QI for proteomics does not perform peptide identifications itself for DDA data. Instead it supports identifications by allowing you to export MS/MS peak lists in formats which can be used to perform peptide searches by various search engines. The resulting identifications can then be imported back into Progenesis QI for proteomics, using a number of different file types, and matched to your detected peptide ions.

**NOTE:** If you have come straight to this stage having performed Automatic Processing then please ensure that you have first filtered your data as described on page 20, to leave only peptide ions with a charge state between 2 and 7. For completeness you may also wish to tag your data as described in the previous sections (pages 33 to 41) before performing the Peptide Identification.

**MS/MS Spectra**

Batch inclusion options for creating export list of ms/ms spectra

Export	Rank	#	Run	Scan number	Exported	Peptide ion intensity	Precursor intensity	(%)	Charge	Precursor m/z	Isotope	Id score	P
<input checked="" type="checkbox"/>	51	2	A1	2686	No	7.6e+007	5.2e+005	0.7	2	539.3007	1		
<input checked="" type="checkbox"/>	33	2	A1	2948	No	7.6e+007	1.5e+006	1.9	2	539.3010	1		
<input checked="" type="checkbox"/>	43	2	A1	2990	No	7.6e+007	8.5e+005	1.1	2	539.3008	1		
<input checked="" type="checkbox"/>	47	2	A1	3040	No	7.6e+007	7.3e+005	1.0	2	539.3007	1		
<input checked="" type="checkbox"/>	58	2	A1	3126	No	7.6e+007	3.0e+005	0.4	2	539.2997	1		
<input checked="" type="checkbox"/>	26	2	A1	2906	No	7.6e+007	2.9e+006	3.8	2	539.3008	1		
<input checked="" type="checkbox"/>	19	2	A1	2855	No	7.6e+007	1.5e+007	20.2	2	539.3004	1		
<input checked="" type="checkbox"/>	44	2	A2	3163	No	7.7e+007	8.2e+005	1.1	2	539.3005	1		
<input checked="" type="checkbox"/>	6	2	A2	2890	No	7.7e+007	6.6e+007	85.5	2	539.3012	1		
<input checked="" type="checkbox"/>	10	2	A2	2931	No	7.7e+007	6.1e+007	79.6	2	539.3006	1		
<input checked="" type="checkbox"/>	17	2	A2	2979	No	7.7e+007	2.5e+007	33.1	2	539.3007	1		
<input checked="" type="checkbox"/>	21	2	A2	3027	No	7.7e+007	7.4e+006	9.7	2	539.3004	1		
<input checked="" type="checkbox"/>	30	2	A2	3077	No	7.7e+007	1.9e+006	2.4	2	539.3007	1		

Peptide ion number 2, m/z 539.3004, retention time 30.518 min, charge +2

Run:A1 Scan number:2686

Retention time (min) vs m/z plot showing peaks at 539, 540, 541, 542.

Intensity vs m/z plot showing a major peak at m/z 539.3004.

Section Complete

**Note:** ensure that **NO** tag filters are applied and that Mascot is the selected search engine.

Determining protein identification is dependent on the availability of MS/MS data for the LC-MS runs. This data may be available but limited if the LC-MS was performed in a data dependant MS/MS detection mode due to under sampling. Under these conditions MS/MS data acquisition is dependent on thresholds and parameters set prior to the acquisition of the LC-MS run.

For this example we are using LC-MS runs containing MS/MS data where the data was acquired in a data dependant mode (DDA).

The Peptide Search page shows the number of MS/MS that have been matched to each peptide ion in the Peptide ion list (see above). MS/MS scans are matched to a peptide ion if their precursor *m/z* and aligned retention time fall within the area of one of the isotopes of the peptide ion. The MS/MS scans which are matched to the displayed peptide ions are shown in the MS/MS spectra list on the right.

The first step is to decide which MS/MS scans you wish to export to be identified. By default this is all the available spectra for the peptide ions displayed in the Peptide ions list. This number is visible on the Export button.


The set can be targeted using the tags and also refined with respect to quantity and quality of the spectra being sent to the search engine.

**Note:** by default the table is ordered on the number of MS/MS spectra available for each peptide ion.

The total number of spectra included in this set is **37276** as shown on the Export button.

Before exporting the spectra, the set can be further refined.

**Note:** many of the abundant peptide ions have a large number of spectra associated with them. To control the number of spectra for each peptide ion, expand the **Batch inclusion options**.



**Batch inclusion options for creating export list of ms/ms spectra**

Rank: **greater than** 5

Peptide ion intensity: less than

Precursor intensity: less than

Precursor intensity (%): less than

Run name: contains

Peptide sequence: contains

Protein accession: contains

Protein description: contains

Peptide ion #: less than

Charge: less than

Scan number: less than

Exported: equal to

Isotope: less than

ID score: less than

Buttons: Include in export, **Exclude from export**, Clear all filters

Export	Rank	#	Run	Scan number	Exported	Peptide ion intensity	Precursor intensity (%)	Charge	Precursor m/z	Isotope	Id score	P
<input checked="" type="checkbox"/>	5	14440	A3	6145	No	1.5e+005	1.0e+005	65.2	3	566.3132	1	
<input checked="" type="checkbox"/>	5	14444	A3	6324	No	1.2e+005	7.4e+004	61.0	2	656.8617	1	
<input type="checkbox"/>	6	1	C3	7200	No	1.1e+008	7.9e+007	71.7	3	805.4420	1	
<input type="checkbox"/>	6	2	A2	2890	No	7.7e+007	6.6e+007	85.5	2	539.3012	1	

Peptide ion number 2, m/z 539.3004, retention time 30.518 min, charge +2

Run: A1 Scan number: 2666

Retention time (min): 29.272, 31.399, 33.529

m/z: 539, 540, 541, 542

Intensity: 0, 5000, 10000, 15000

m/z: 200, 400, 600, 800, 1000

Section Complete

**For example:** We will make use of the 'Rank' value to reduce the number of Spectra being used for each peptide ion to a maximum of 5.

The 'Rank' of each MS/MS spectra is determined by comparing its % value against all other spectra matched to the same peptide ion.

Export	Rank	#	Run	Scan number	Exported	Peptide ion intensity	Precursor intensity (%)	Charge	Precursor m/z	Isotope	Id s
<input checked="" type="checkbox"/>	5	801	C2	3577	No	5.3e+005	1.8e+005	33.7	2	684.8735	1
<input checked="" type="checkbox"/>	5	404	A3	2362	No	2.9e+005	3.4e+005	13.9	3	337.0436	1
<input type="checkbox"/>	6	1800	A3	4618	No	2.9e+005	1.1e+005	39.4	2	783.8677	1
<input type="checkbox"/>	6	1756	C2	8216	No	2.9e+005	1.6e+005	54.4	2	729.8898	1

The rank of each MS/MS spectrum found by comparing its '%' values against all other spectra matched to the same peptide ion.

**Note:** the % value for each spectra is the Precursor intensity as a percentage of the Peptide ion intensity

Set the Rank filter to 'greater than' 5 and click **Exclude from export** this reduces the number to spectra to export to **18131**

Limiting the 'fragment ion count' (FIC) for the spectra being exported can improve the quality of the spectral data being used in the search by removing noisy peaks.

For example for the current spectra, reduce the FIC to 40.

**Note:** the effect this has on the number of peaks in the spectra. This 'limitation' is applied to all the spectra being exported; hence the export file size will be reduced.

The image displays two screenshots of the Progenesis QI software interface, illustrating the effect of limiting the fragment ion count (FIC) on the resulting MS/MS spectra.

**Top Screenshot:** The 'MSMS Preprocessing' section shows the 'Limit fragment ion count' checkbox is **unchecked**. The resulting spectrum (Run: C2 Scan number: 3762) shows a high number of peaks, indicating a high FIC.

**Bottom Screenshot:** The 'MSMS Preprocessing' section shows the 'Limit fragment ion count' checkbox is **checked**. The resulting spectrum (Run: C2 Scan number: 3762) shows a significantly reduced number of peaks, indicating a lower FIC. A large black arrow points from the top spectrum to the bottom spectrum, highlighting the reduction in peak count.

For this example, we will **NOT** limit the fragment count, so leave it un-ticked (the default setting).

## Performing an MS/MS Ion Search

Having chosen 18131 spectra to export, as described above:


1. Select appropriate search engine i.e. Mascot
2. Click 'Export 18131 ms/ms spectra' to save search as file
3. Perform search on appropriate search engine and save results file
4. Click 'Import search results', locate results file and open

Please refer to **Appendix 12A (page 98)** for details of the 'Search Engine' parameters used in this example

**Tip:** For exporting mzIdentML of results to PRIDE where the Search Method is MascotSearch you **must select** additional **Mascot Export Fields** when exporting your search results from Mascot to Progenesis: (Protein sequence (Protein Hit Information) and Start and End (Peptide Match Information)), **Appendix 12A (page 98)**

**Note:** the blue link provides you with details on the appropriate formats for exporting search results and access to additional formats

**Note:** an example Search Results file, from a MS/MS Ion search, is available in the folder you restored the Archive to (Search Results.xml). Select the 'Mascot' method and import this file to see results like those below.



The image shows a screenshot of the Progenesis QI software interface during the 'Identify Peptides' workflow step. The 'Identify Peptides' tab is active, and a red box highlights the 'Identify Peptides' button in the top navigation bar. The main window displays a table of identified peptides with columns for '#', 'MS/MS', 'Proteins', 'Tag', and 'Score'. A red box highlights the first few rows of the table. Below the table, there are options for 'Select the search program you're using:' (Mascot) and 'Export 0 ms/ms spectra'. The 'MSMS Preprocessing' section is also visible, with a checkbox for 'Limit fragment ion count' set to 40 and a checkbox for 'Deisotoping and charge deconvolution'. A red box highlights the 'Import search results' button. A dialog box titled 'QIP Import search results' is open, displaying the message: '8898 spectrum matches have been imported and assigned to peptide ions'. Below the dialog, there is a table of peptide ions with columns for 'Peptide ion number', 'm/z', 'Retention time', 'Charge', 'Precursor m/z', 'Isotope', 'Id score', and 'P'. A red box highlights the first row of this table. At the bottom, there is a plot showing 'Retention time (min)' vs 'm/z' for 'Run: C2 Scan number: 3762'. The plot shows a peak at m/z 901.462 with a retention time of 34.662 min. A red box highlights the peak. The 'Section Complete' indicator is visible at the bottom right.


On importing the Search results the Peptide ions table updates to reflect the identified proteins and the relevant score for each searched peptide ion.

In order to review the quality of the **Peptide Search** results click on the next stage in the workflow, **QC Metrics**.

## Stage 10: QC Metrics

Progenesis QI for proteomics includes a number of batch-based QC metrics which you can examine to increase confidence in, or identify issues with, your data. These views are presented at the “QC Metrics” page. By setting up experimental groupings that correspond to your batch metadata, you are able to investigate whether there are any systematic effects associated with your processing. This gives you confidence that your conclusions are not affected by technical biases within the course of the experiment.

The charts are updated with changes in your active experimental design, allowing you to examine the effects of all stages of your processing if you wish through using different batch divisions.



### Metrics

There are two sets of metrics. Firstly, the overview page provides nine at-a-glance summary measurements which cover the experiment as a whole. Then, the remainder of the pages provide more detailed information and bring in the batch-by-batch detail, some of which builds on the overview charts:

- Sample preparation metrics – highlighting issues or problems with the preparation of your samples: Missed cleavages, modifications and abundance dynamic range
- Instrument metrics – highlighting whether your chromatography column and mass spectrometer are configured and performing correctly: Mass accuracy, abundance dynamic range (again), precursor charges, MS1 scan rates
- Experiment metrics – concerning the identified proteins and peptides in your experiment, allowing you to pick out any outlying runs or conditions: Proteins, peptides, peptides per protein, % of peptide ions identified, proteins per condition

More detailed information on the QC metrics is available on the FAQ pages.


### Interpretation and use

It's worth noting that there is no simple 'right' or 'wrong' answer as to whether your data are of high enough quality.

**Note:** QC measurements are designed to flag up potential issues for investigation, hence there must be an assignable cause to the variation observed to render QC practical.

Furthermore, the QC metrics will be rebuilt if you delete and re-do any stages leading up to them (for example, clearing identifications and re-searching with changed identification parameters, or altering the alignment). Hence, if you archive your experiment before and after your changes, or export the charts before and after, you can see the effects of your changes upon the quality measurements – a very useful method for assessing your interventions.

**Note:** that the overview metrics can each be expanded by clicking on the top-right icon in the sub-window. Additionally, hovering over a column will bring up a tooltip containing quantitative information on the results.



You can also tag metrics of interest or concern with a comment for your records, which is saved with the experiment. To do this, click the empty flag icon in the tab header for the given metric (or in the tile header on the overview screen):



This will pop up a dialog allowing you to enter a comment describing why this metric has been flagged.

## Reporting

You can export a report for your own records, or to pass on to another member of your team:

### Report all metrics:

Generates a report containing all metrics shown in the application (overview, sample preparation metrics, instrument metrics and experiment metrics). This may be useful for documentation purposes, or to verify the quality of your experiment.

### Report flagged metrics:

Generates a report of only metrics you have flagged, along with the message you provided. This may be useful for giving to a technician or other team member, to highlight areas of the experiment that need improvement.



In order to review and refine the quality of the **Peptide Search** results click on the next stage in the workflow, **Refine Identifications**.

## Stage 11: Refine Identifications

In this example the organism under study is *Clostridium difficile*

**Note:** before removing any identifications, make sure there are **no** tag filters applied at the Identify peptides stage.

As an **example** the following section describes how sequential filtering of the Peptide results can be performed using the following thresholds described below:

- Remove identifications with a Score less than 40
- Remove identifications where less than 2 hits were returned
- Remove all identifications where the Protein Description **Contains** 'hypothetical'
- Remove all identifications where the Protein Description **Contains** the following: 'Like', 'Putative', 'Probable', 'Potential', 'Predicted' and 'Partial'
- Remove all identifications where the Protein Description **Doesn't contain** 'Clostridium difficile'


The screenshot shows the 'Refine Identifications' window in Progenesis QI. The 'Batch deletion criteria' panel is active, with 'Score' set to 'less than 40'. The 'Delete matching search results' checkbox is checked. The main table displays a list of peptide identifications with columns for #, Score, Hits, m/z, RT (mins), Charge, Mass, Mass err, Sequence, Accession, Modifications, and Description. The bottom of the table shows '2903 search results. 512 matching batch delete options.' A 'Section Complete' button is visible in the bottom right corner.

On the Batch detection options panel, set the Score to less than 40, then **Delete matching search results**.

**Note:** the search results matching the filter criteria turn pink and the total is displayed at the bottom of the table (512 matching out of 2903)

**Note:** a dialog warns you of what you are about to delete

Click **Yes**.



Now click **Reset the criteria** to clear the previous threshold then apply the next filter 'Hits: less than 2'.

Then in the 'Description' field enter 'hypothetical' and delete matching search results. Then also in Description field enter the 'regular expression': **regex: like|puta|prob|pote|pred|part** and delete matching search results.

**Refine Identifications**

If your peptide identifications include unwanted or irrelevant results, you can remove them here.

- Specify a set of deletion criteria**  
In the batch deletion criteria, enter the property values for a set of identifications you want to delete.
- Delete the unwanted identifications**  
To delete the identifications you don't want, click either:
  - Delete Matching Search Results, to delete the highlighted IDs
  - Delete Non-matching Search Results, to delete the IDs that are not highlighted
- Reset the criteria to start again**  
To specify another batch of identifications to delete, click Reset the Criteria and then return to step 1 above.

Batch deletion criteria:

- Score: less than
- Hits: less than
- Mass: less than
- Absolute mass error (ppm): less than
- m/z: less than
- Retention Time: less than
- Sequence Length: less than
- Charge: less than
- Sequence: contains
- Accession: contains
- Description: contains **regex: like|puta|prob|pote|pred|part**
- Modifications: contains

Delete matching search results  Delete non-matching search results

**Delete 714 search results?**  
Are you sure you want to permanently delete 714 peptide search results?

#	Total Hits	m/z	RT (min)	Charge	Tag
1	10	805.4410	55.59	3	
2	60	539.3004	30.52	2	
3	55	598.3214	31.49	2	
4	10	624.2913	28.46	2	
5	35	604.3377	27.34	2	
6	15	626.3146	39.56	2	
7	10	462.2705	35.00	2	
8	12	753.8284	31.44	2	
9	30	656.8613	44.60	2	
10	18	663.8693	47.56	2	
11	20	595.3191	37.47	2	
12	18	573.8033	25.39	2	
14	60	573.3242	42.39	2	

2065 search results. 714 matching batch delete options.

Finally alter the Description to 'doesn't contain Clostridium difficile' and delete the matching search results.

**Refine Identifications**

If your peptide identifications include unwanted or irrelevant results, you can remove them here.

- Specify a set of deletion criteria**  
In the batch deletion criteria, enter the property values for a set of identifications you want to delete.
- Delete the unwanted identifications**  
To delete the identifications you don't want, click either:
  - Delete Matching Search Results, to delete the highlighted IDs
  - Delete Non-matching Search Results, to delete the IDs that are not highlighted
- Reset the criteria to start again**  
To specify another batch of identifications to delete, click Reset the Criteria and then return to step 1 above.

Batch deletion criteria:

- Score: less than
- Hits: less than
- Mass: less than
- Absolute mass error (ppm): less than
- m/z: less than
- Retention Time: less than
- Sequence Length: less than
- Charge: less than
- Sequence: contains
- Accession: contains
- Description: **doesn't contain** **clostridium difficile**
- Modifications: contains

Delete matching search results  Delete non-matching search results

**Delete 177 search results?**  
Are you sure you want to permanently delete 177 peptide search results?

#	Total Hits	m/z	RT (min)	Charge	Tag
1	10	805.4410	55.59	3	
2	60	539.3004	30.52	2	
3	55	598.3214	31.49	2	
4	10	624.2913	28.46	2	
5	35	604.3377	27.34	2	
6	15	626.3146	39.56	2	
7	10	462.2705	35.00	2	
8	12	753.8284	31.44	2	
9	30	656.8613	44.60	2	
10	18	663.8693	47.56	2	
11	20	595.3191	37.47	2	
12	18	573.8033	25.39	2	
14	60	573.3242	42.39	2	

1351 search results. 177 matching batch delete options.

Having applied all the filters there will be **1174** search results remaining

To validate the Peptide search results at the protein level click on **Resolve Conflicts** (bottom left).



## Resolving Conflicts

This stage allows you to examine the behaviour of the identified peptides and choose to resolve any conflicts for the various peptide assignments at the protein level.

The **Resolve Conflicts** stage is now accessed from the bottom left of the **Refine Identifications** stage.


The number of conflicts you have to resolve will depend on the scope and stringency of the filters you apply at the **Refine Identifications** stage.

Details on performing conflict resolution are available in Appendix 8 (page 87)

**For this guide: we will NOT resolve the conflicts.**

The following pages describe the handling of **Protein grouping** and **Protein Quantitation**. Options for these are available at both **Resolve Conflicts** and **Review Proteins**

**Note:** the default settings are to **Use Protein Grouping** and **Relative Quantitation using HiN** (peptides =3)



## Protein Grouping

Where proteins are identified containing the same peptides then they are effectively indistinguishable aside from the score. Also where one protein contains only peptides that represent a subset of another protein's peptides, the protein with fewer peptides can be subsumed into a group where the protein with the greater number is deemed the 'patriarch' and appears in the proteins table.

With protein grouping switched on (default setting) protein groups, the additional members are indicated by a bracketed number located after the Accession number. Taking **flagellin subunit** as an example, there are 9 additional group members. There are no conflicts and all the peptides are uniquely assigned.

Accession	Peptides	Unique	Conflict	Score
gj 260682215 (+1)	24	11	18	2.081
gj 126700407	16	16	0	1.641
gj 254973900 (+9)	15	15	0	1.951
gj 126698450 (+1)	15	15	0	1.951
gj 126697969 (+1)	10	10	0	9.2

#	Σ	Score	Hits	Correlation	Mass	Mass error (p...)	RT (mins)	Charge	Tag
67	2	63.9	5	0.965	1669.889	-0.0445	42.9	3	
121	1	103	5	0.963	1669.888	-0.562	42.9	2	
148	2	101	5	0.921	1230.609	-0.423	24	2	
172	2	60.9	3	0.951	2317.115	0.2	39.5	3	

Protein: gj|254973900 flagellin subunit [Clostridium difficile QCD-66c26]

No protein selected

Peptide ion views Protein resolution

Conflicting proteins

Accession	Peptides	Unique	Conflict	Protein Score
-----------	----------	--------	----------	---------------

Peptide ions of conflicting protein

#	Σ	Score	Hits	Correlation	Mass	Mass error (p...)	RT (mins)	Charge	Tag	Abu
---	---	-------	------	-------------	------	-------------------	-----------	--------	-----	-----

Refine Identifications Protein options... Section Complete

Now open the protein options and untick **Use protein grouping** to examine the effect of the process.

As grouping is switched off the grouped proteins appear with conflicts to the other group members and there will be no unique peptides for each group member.

The screenshot shows the Progenesis QI software interface. At the top, the 'Experiment design' is set to 'AC'. Below this, there are two main panels: 'Proteins' and 'Peptide ions of gj|260682017'. The 'Proteins' panel shows a list of proteins with columns for Accession, Peptides, Unique, Conflict, Score, and Tag. A red box highlights a group of proteins where the 'Unique' column is 0 for all entries. Below the protein list, there are two protein descriptions: 'Protein: gj|260682017 flagellin subunit [Clostridium difficile CD196]' and 'Protein: gj|254973900 flagellin subunit [Clostridium difficile QCD-66c26]'. The 'Peptide ions of gj|260682017' panel shows a list of peptide ions with columns for #, Σ, Score, Hits, Correlation, Mass, Mass error (p...), RT (mins), Charge, Tag, Abundance, and Conf. Below this, there are two more panels: 'Conflicting proteins for peptide ion 452' and 'Peptide ions of gj|254973900'. The 'Conflicting proteins for peptide ion 452' panel shows a list of proteins with columns for Accession, Peptides, Unique, Conflict, Protein Score, and Tag. A red box highlights a group of proteins where the 'Unique' column is 0 for all entries. The 'Peptide ions of gj|254973900' panel shows a list of peptide ions with columns for #, Σ, Score, Hits, Correlation, Mass, Mass error (p...), RT (mins), Charge, Tag, Abundance, and Cor. At the bottom of the interface, there are buttons for 'Refine Identifications', 'Protein options...', and 'Section Complete'.

All the group members will appear in both tables with large numbers of conflicts with the other group members.

**Note:** flagellin has **no unique** peptides as they are all present in flagellin subunit protein hence the reason for grouping. As a result, all the conflicts are internal to the group.

Now set the **Protein Options** back to **Use protein grouping**.

The screenshot shows the 'QIP Protein quantitation options' dialog box. The 'Quantitation method' is set to 'Relative Quantitation using Hi-N'. The 'Number of peptides to measure per protein (N):' is set to 3. The 'Use protein grouping' checkbox is checked, and the text next to it reads 'i.e. hide proteins whose peptides are a subset of another protein's.'. At the bottom of the dialog box, there are 'OK' and 'Cancel' buttons.

## Protein Quantitation options

There are 4 options with which to control how the Protein Quantification is performed by Progenesis QI for proteomics.

The default option that will be applied is **Relative Quantitation using Hi-N**.

If you have selected one of the other options during the setup of the Auto Processing of your data (Stage 1) then this option will be applied.





Hi-N in Progenesis QI for proteomics is an implementation of Hi-3 as described by Silva *et al.* [References]. After peptide and protein identification, the abundance of each peptide is calculated from all its constituent peptide ions.

For each protein, the N most abundant peptides (N being set according to the user selection) have their **abundances averaged** to provide a reading for the protein signal.

The ranking of peptide abundance is based on the integrated value across all the runs, allowed by the accurate alignment and lack of missing values. This gives added confidence in the peptide selection, taking all runs into account to make the ranking robust.

The averaged abundance readings not only make possible the **relative quantitation** of the same protein across all runs to be determined but also allow, with the inclusion of a known amount of a calibrant protein in each run, this to be converted to an **absolute** reading for protein amount.

To generate values for **absolute quantitation** enter the accession number (i.e. for enolase) and amount for the calibrant.



The absolute amounts, based on the calibrant used, are reported at the Review Proteins stage as additional columns (one for each condition) following the protein description.

**Review Proteins**  
Using this screen, you can find the proteins of interest in your experiment.

No filter applied

Max fold change	Highest Mean	Lowest Mean	Description	Amount (fmol) - A	Amount (fmol) - C
8.6	A	C	transketolase, central and C-terminal (Sedoheptulose-7-phosphate-D-glyceraldehyde-3-phosphate glycolaldehy	20.1	2.34
3.54	C	A	fructose-1,6-bisphosphate aldolase [Clostridium difficile 630]	2.13	7.54
11.1	A	C	thioredoxin 2 (Trx2) [Clostridium difficile 630]	482	43.3
4.97	A	C	transcription elongation factor GreA [Clostridium difficile 630]	90.1	18.1

The ranking of peptide abundance is based on the integrated value across all the runs, allowed by the accurate alignment and lack of missing values. This gives added confidence in the peptide selection, taking all runs into account to make the ranking robust. These Hi-N methods allow the relative and absolute comparison of proteins within the same run.

**Note:** When there are **peptide conflicts** (a peptide is shared between two proteins, for example) it is important to assign the signal correctly for absolute quantitation. To do this, Progenesis QI for proteomics carries out a two-step process. Firstly, Hi-N is carried out only on the N most abundant **unique** (non-conflicting) peptides of the proteins concerned. This provides a ratio estimate for the two proteins based only on unshared peptides. The abundance of any shared peptides is then divided and allotted in this ratio between the two proteins, and the full Hi-N calculation is then applied using the divided values for conflicted peptides.

Naturally, if the conflicted peptides are not among the N most abundant in either protein initially, then this will not cause any difference in the result.

To generate values for **Absolute Quantitation for HCP** (Host Cell Proteins) select the 5<sup>th</sup> option on the drop down at the beginning of this section.

**Note:** there are 2 differences between this HCP-specific option and the Absolute Quantitation option:

- The mass of protein present is also reported in ng in addition to the amount in fmol
- The calculation of protein amounts and masses present is based on measurements of the calibrant protein present in the same run, and not a pooled measurement over all the runs. It does still assume a specified amount of a calibrant in every run. However, it does not pool information across the runs to derive its relationship between the calibrant and observed abundance. The amount of a contaminant is instead calculated directly by relation to the observed amount of the calibrant in the same run.

**Review Proteins**  
Using this screen, you can find the proteins of interest in your experiment.

1 **Set the quantitation options**  
If you've not already done so, choose between relative and absolute quantitation, use of Hi-N, protein grouping and more.

No filter applied    Create...    Search

Description	Amount (fmol) - A	Amount (fmol) - C	Amount (ng) - A	Amount (ng) - C
formate--tetrahydrofolate ligase [Clostridium difficile 630]	9.45	11.2	0.57	0.674
fructose-1,6-bisphosphate aldolase [Clostridium difficile 630]	4.22	5.07	0.14	0.168
gamma-aminobutyrate metabolism dehydratase/isomerase [Clostridium difficile 630]	31.9	14.4	1.78	0.802
glyceraldehyde-3-phosphate dehydrogenase [Clostridium difficile 630]	6.56	14.4	0.238	0.522

As for **Absolute Quantitation** enter the accession number and amount for the calibrant.

**QIP Protein quantitation options**

Quantitation method:  
Absolute Quantitation for HCP using Hi-N

Number of peptides to measure per protein (N): 3

Calibrant accession:   
Amount (fmol):

Use protein grouping i.e. hide proteins whose peptides are a subset of another protein's.

OK Cancel

**QIP Protein quantitation options**

Quantitation method:  
Absolute Quantitation for HCP using Hi-N

Number of peptides to measure per protein (N): 3

Calibrant accession: **gij126700790**  
Amount (fmol): **50**

Use protein grouping i.e. hide proteins whose peptides are a subset of another protein's.

OK Cancel

The absolute amounts, based on the calibrant used, are reported at the Review Proteins stage as additional columns (one for each condition) following the protein description.


**How the values for the calibrant behave across runs comparing Absolute Quant and Absolute Quant for HCP using the same calibrant**

Description	Amount (fmol) - A	Amount (fmol) - C
elongation factor Ts [Clostridium difficile 630]	46.5	5.06
enolase (2-phosphoglycerate dehydratase) (2-phospho-D-glycerate hydro-lyase) [Clostridium difficile 630]	25.4	74.6
FOF1 ATP synthase subunit beta [Clostridium difficile 630]	3.31	11.4


Description	Amount (fmol) - A	Amount (fmol) - C	Amount (ng) - A	Amount (ng) - C
elongation factor Ts [Clostridium difficile 630]	91.9	3.4	3.06	0.113
enolase (2-phosphoglycerate dehydratase) (2-phospho-D-glycerate hydro-lyase) [Clostridium difficile 630]	50	50	2.31	2.31
FOF1 ATP synthase subunit beta [Clostridium difficile 630]	6.56	7.67	0.327	0.382

**Relative Quantitation** can also be performed comparing a single protein across all the runs using only the unique or non-conflicting peptides. Select the third method from the drop down

Using non conflicting peptides



The relative Quantitation can also be performed using all peptides.



**Note:** if you have performed conflict resolution then there will be no difference between these methods.

Now move to the Review Proteins section by clicking on **Review Proteins** icon on the workflow at the top of the screen.

## Stage 12: Review Proteins

The **Review Proteins** stage opens displaying details for all the proteins. You can now create tags at the protein level for proteins with an Anova p value  $\leq 0.05$  and Max Fold change  $\geq 2$  by right clicking on the table and selecting **Quick Tags**...

**Review Proteins**

Using this screen, you can find the proteins of interest in your experiment.

1 **Set the quantitation options**  
If you've not already done so, choose between relative and absolute quantitation, use of Hi-N, protein grouping and more.

2 **Create a shortlist to review**  
In the table, sort and filter the proteins based on their measurements, to generate a shortlist for further review.

3 **Review the proteins**  
For each protein of interest, review its peptide measurements and correlations.

4 **Export data for further processing**  
By exporting your data to external tools, there's no limit to your analysis.

Accession	Peptides	Unique peptides	Confidence score	Anova (p)	q Value	Tag	Max fold change	Highest Mean	Lowest Mean	Description
gi 126699603	3	3	236	1.23E-06	3.26E-06		5.37	A	C	decarboxylase
gi 126701103	1	1	51.7	5.06E-06	1E-05		6.17	A	C	ribose-5-phospho- transketolase, c
gi 126699939	3	3	236	5.3E-06	1E-05		8.6	A	C	fructose-1,6-bi
gi 126697972	2	2	130	6.61E-06	1.09E-05		3.54	A	C	thioredoxin 2
gi 126699971	5	5	296	1.33E-05	1.33E-05		11.1	A	C	transcription el
gi 126701179	3	3	236	1.33E-05	1.33E-05		4.97	A	C	cell surface pro
gi 254974053	1	1	73.3	1.41E-05	1.41E-05		898	A	C	elongation fact
gi 126699756	3	3	236	1.41E-05	1.41E-05		898	A	C	hemagglutinin
gi 260682215 (+1)	24	11	73.3	1.41E-05	1.41E-05		898	A	C	flagellin subun
gi 254973900 (+9)	15	15	73.3	1.41E-05	1.41E-05		898	A	C	ferredoxin/Flav
gi 126697690	6	6	298	2.99E-05	2.99E-05		7.16	A	C	PTS system ma
gi 126700634	3	3	73.3	3.46E-05	2.09E-05		7.16	A	C	transketolase, f
gi 126699940	1	1	73.3	3.46E-05	2.09E-05		7.16	A	C	transketolase, f
gi 126700297	2	2	164	3.67E-05	2.7E-05		2.29	A	C	propanediol ut
gi 126697583	1	1	79.3	4.33E-05	3E-05		3.78	A	C	DNA binding p
gi 126700790 (-1)	8	8	744	4.75E-05	3E-05		2.94	A	C	enolase 12-pho

**Selected protein: thioredoxin 2 (Trx2) [Clostridium difficile 630]**

View peptide measurements

Standardised Expression Profiles

Peptide profiles: ■ Contributes to protein measurements ■ Does not contribute to protein measurements

As an example, let us explore Thioredoxin 2. The table indicates that this protein is most highly expressed in Condition A by 11.1 fold over the lowest condition (C).

To view the corresponding peptide measurements for the current protein either double click on the protein in the table or use the **View peptide measurements** below the table.

**Review selected protein**

Review the selected protein's identified peptides and validate their expression patterns.

1 **Choose the level of detail**  
View the properties and expression profiles of either peptides or individual peptide ions:

2 **Compare expression profiles**  
Select peptides in the table to show their expression profiles in the chart below.

3 **Resolve any quantitative outliers**  
Tag any peptide ions whose expression profile is an outlier for this protein.

**Experiment design**  
Review your data from a different perspective:

Current design: AC

Quantifiable proteins displayed: 145

**Review selected protein**

Accession: gi|126699971  
Description: thioredoxin 2 (Trx2) [Clostridium difficile 630]

Identifier	Ions	Score	Correlation	Anova (p)	Max Fold Change	Highest Mean	Lowest Mean	Tag	Abundance	Neutral Mass
48.79_1204.6891n	1	69.090	0.976	0.000188	15.6	A	C		2.266E+06	1204.6891
42.41_2130.0509n	2	66.320	0.995	5E-06	8.3	A	C		1.948E+06	2130.0509
40.48_1220.6834n	1	48.680	0.998	7.52E-05	11.5	A	C		4.364E+05	1220.6834
46.97_1458.7017n	1	60.280	0.999	0.000972	16.8	A	C		2.27E+05	1458.7017
33.27_2146.0462n	2	44.060	0.985	6.33E-05	13	A	C		2.784E+05	2146.0462

Standardised Expression Profiles

Peptide profiles: ■ Contributes to protein measurements ■ Does not contribute to protein measurements

The solid icon in  $\Sigma$  (click to order on) column indicates that the peptide contributes to protein measurements. **Note:** a correlation score for each peptide is available to facilitate the validation of peptide expression.

By selecting all the peptides you can compare the pattern of expression across all the samples allowing you to identify any 'atypical' behaviour assigned to the current protein.

**Review selected protein**  
 Review the selected protein's identified peptides and validate their expression patterns.

1 **Choose the level of detail**  
 View the properties and expression profiles of either peptides or individual peptide ions:  
 Show:  Peptides  Peptide ions  
 Tip: you can also double-click a peptide to select and view its component ions.

2 **Compare expression profiles**  
 Select peptides in the table to show their expression profiles in the chart below.

3 **Resolve any quantitative outliers**  
 Tag any peptide ions whose expression profile is an outlier for this protein.

You can then review their identifications in more detail at the [Refine Identifications](#) step.

**Experiment design**  
 Review your data from a different perspective:  
 Current design:

**Accession:** gil126699971  
**Description:** thioredoxin 2 (Trx2) [Clostridium difficile 630]

Fold Change	Highest Mean	Lowest Mean	Tag	Abundance	Neutral Mass	Retention Time (mins)	Peptide Sequence	Modifications
A	C	C		2.266E+06	1204.6891	48.79	VLGLPTMAIYK	
A	C	C		1.948E+06	2130.0509	42.409	VDEVTKDDATVPNIENMIK	
A	C	C		4.364E+05	1220.6834	40.481	VLGLPTMAIYK	[7] Oxidation (M)
A	C	C		2.27E+05	1458.7017	46.97	DDATVPNIENMIK	
A	C	C		2.784E+05	2146.0462	33.269	VDEVTKDDATVPNIENMIK	[17] Oxidation (M)

**Standardised Expression Profiles**

Peptide profiles:  Contributes to protein measurements  Does not contribute to protein measurements

Scroll to the right on the table for the sequence and modification details. Use the navigation buttons on the top left either to review each protein's peptides by stepping through the list or return to the protein list and select a specific protein to review.

**Note:** as the quantitation is based on the top 3 peptides, the remainder are shown grey on the graph

The behaviour of individual peptide ions can also be reviewed by selecting **peptide ions** as the level of detail to review.

**Review selected protein**  
 Review the selected protein's identified peptides and validate their expression patterns.

1 **Choose the level of detail**  
 View the properties and expression profiles of either peptides or individual peptide ions:  
 Show:  Peptides  Peptide ions  
 Tip: you can also double-click a peptide to select and view its component ions.

2 **Compare expression profiles**  
 Select peptide ions in the table to show their expression profiles in the chart below.

3 **Resolve any quantitative outliers**  
 Tag any peptide ions whose expression profile is an outlier for this protein.

You can then review their identifications in more detail at the [Refine Identifications](#) step.

**Experiment design**  
 Review your data from a different perspective:  
 Current design:

**Accession:** gil126699971  
**Description:** thioredoxin 2 (Trx2) [Clostridium difficile 630]

an	Lowest Mean	Tag	Abundance	m/z	Charge	Retention Time (mins)	Mass error (ppm)	Peptide Sequence	Modifications
C			2.27E+05	730.3581	2	46.97	-0.51	DDATVPNIENMIK	
C			1.384E+06	711.0243	3	42.429	-0.17	VDEVTKDDATVPNIENMIK	
C			5.637E+05	1066.0326	2	42.39	-0.34	VDEVTKDDATVPNIENMIK	
C			2.14E+05	716.3563	3	33.269	0.25	VDEVTKDDATVPNIENMIK	[17] Oxidation (M)
C			6.443E+04	1074.03	2	33.269	-0.48	VDEVTKDDATVPNIENMIK	[17] Oxidation (M)
C			2.266E+06	603.3518	2	48.79	0.08	VLGLPTMAIYK	
C			4.364E+05	611.349	2	40.481	-0.42	VLGLPTMAIYK	[7] Oxidation (M)

**Standardised Expression Profiles**

Peptide ion profiles:  Contributes to protein measurements  Does not contribute to protein measurements

**Note:** peptide ions with different charge states are now displayed, these are deconvoluted at the peptide level.

Returning to the protein level, with protein grouping switched on (default setting) the additional members are indicated by a bracketed number located after the Accession number. The additional members can be seen by holding the cursor over the accession number.

**Review Proteins**

Using this screen, you can find the proteins of interest in your experiment.

- Set the quantitation options**  
If you've not already done so, choose between relative and absolute quantitation, use of Hi-N, protein grouping and more.
- Create a shortlist to review**  
In the table, sort and filter the proteins based on their measurements, to generate a shortlist for further review.
- Review the proteins**  
For each protein of interest, review its peptide measurements and correlations:
- Export data for further processing**  
By exporting your data to external tools, there's no limit to your analysis.

Accession	Peptides	Unique peptides	Confidence score	Anova (p)	q Value	Tag	Max fold change	Highest Mean	Lowest Mean	Description
gi 126700198	1	1	46.6	0.000104	4.9E-05		9.14	A	C	phosphatase, 2C
gi 126700129	3	3	267	0.000125	5.68E-05		2.32	A	C	translation inhibi
gi 254974605	4	4	311	0.000154	6.79E-05		2.78	C	A	electron transfer
gi 254973854 (+4)	5	5	428	0.000166	7.08E-05		3.03	C	A	60 kDa chaperon
gi 126698915								C	A	30S ribosomal pr
gi 126698631 (+)								C	A	cell wall-binding
gi 255101963 (+)								A	C	cell surface prote
gi 126699972								A	C	thioredoxin redu
gi 126699583	2	2	197	0.000219	8.03E-05		5.1	C	A	RNA-binding pro
gi 126697654	3	3	215	0.00023	8.22E-05		2.45	A	C	30S ribosomal pr
gi 126699299	2	2	311	0.000248	8.44E-05		3.34	A	C	dinitrogenase iro

**Selected protein: 60 kDa chaperonin [Clostridium difficile QCD-66c26]**

View peptide measurements

Section Complete

To view members peptides click on **View peptide measurements**.

Click on **Select peptides of** to show the list of additional group members. Then as you click on one of the member proteins its peptides are highlighted in the table.

**Review selected protein**

Review the selected protein's identified peptides and validate their expression patterns.

Accession: gi|254973854

Description: 60 kDa chaperonin [Clostridium difficile QCD-66c26]


#	Score	Correlation	Anova (p)	Max Fold Change	Highest Mean	Lowest Mean	Tag	Abundance	m/z	Charge	Retention Time (mins)
770	66.390	0.971	0.00273	2.07	C	A		1.784E+05	619.8642	2	32.549
862	58.830	0.982	3.81E-05	3.52	C	A		2.952E+05	766.9192	2	40.786
978	77.730	0.988	0.000371	2.99	C	A		2.918E+05	696.7028	3	54.573
2005	106.050	0.991	6.88E-05	3.83	C	A		2.153E+05	1044.5499	2	54.624
1185	68.460	0.997	1.17E-06	51.3	C	A		1.44E+05	535.8258	2	35.849
13282	50.250	0.755	0.0187	10.6	C	A		4074	717.079	3	72.049

**Standardised Expression Profiles**


Section Complete




You can tag all the peptide ions for a protein or edit and tag accordingly for 'atypical behaviour'. Then remove these peptides at the **Refine identification** stage in the Workflow. Example, at the Protein level of **Review proteins** use the search facility to locate the protein with accession gj126700790.



Then clear the search box (click on 'x') leaving the protein highlighted in the table. Double click on the protein and then click on Show: **Peptide ions** then select 2 of the poorly correlated right click on the selection and Create a New tag, **Poorly Correlated**.



To delete the tagged peptide ions (Poorly Correlated). Click on **Refine Identifications**, then set up a tag filter to restrict the left hand table to only those peptide ions tagged with Poorly Correlated.



Then highlight the two entries in the left-hand table, the corresponding rows will be selected on the right-hand table. Right click on the highlighted rows and Delete selected peptide(s)

**Refine Identifications**

If your peptide identifications include unwanted or irrelevant results, you can remove them here.

- Specify a set of deletion criteria**  
In the batch deletion criteria, enter the property values for a set of identifications you want to delete.
- Delete the unwanted identifications**  
To delete the identifications you don't want, click either:
  - Delete Matching Search Results, to delete the highlighted IDs
  - Delete Non-matching Search Results, to delete the IDs that are not highlighted
- Reset the criteria to start again**  
To specify another batch of identifications to delete, click **Reset the Criteria** and then return to step 1 above.

**Batch deletion criteria**

Score: less than  
Hits: less than  
Mass: less than  
Absolute mass error (ppm): less than  
m/z: less than  
Retention Time: less than

Sequence Length: less than  
Charge: less than  
Sequence: contains  
Accession: contains  
Description: contains  
Modifications: contains

Delete matching search results  Delete non-matching search results

#	Score	Hits	m/z	RT(mins)	Charge	Mass	Mass err.	Sequence	Accession	Modifications
4492	---	---	652.97	41.13	3	1955.8	-0.07	NSDLNTVS...	gi 1267004...	cell wall binding protein [Clostridium
4521	---	---	600.30	47.32	3	1797.8	1.12	FIEEIGYVN...	gi 1266988...	30S ribosomal protein S16 [Clostrid
4593	53.44	4	700.87	18.38	2	1399.7	-0.70	IAEELGNRE...	gi 2549763...	cell surface protein [Clostridium dif
4593	53.44	4	700.87	18.38	2	1399.7	-0.70	IAEELGNRE...	gi 3575510...	Erk/YbiS/YctS/YnhG [Clostridium d
4625	105.85	5	974.01	61.34	2	1946.0	-0.04	VNQGITTE...	gi 2964520...	enolase [Clostridium difficile NAP0
4625	105.85	5	974.01	61.34	2	1946.0	-0.04	VNQGITTE...	gi 1267007...	enolase (2-phosphoglycerate dehy
4634	---	---	1116.5	30.02	2	2			1267001...	translation inhibitor endoribonucle
4642	60.64	5	863.41	22.56	2	1724.8	-0.26	QADREGYP...	gi 1266991...	rubrythrin [Clostridium difficile 63
4679	58.76	3	540.29	62.31	3	1617.8	0.34	DAQIAEVV...	gi 1267006...	PTS system mannose-specific transp

**Tag filter applied**  
peptide ions may be hidden

#	Total Hits	m/z	RT (min)	Charge	Tag
4625	10	974.001	61.34	2	
10041	---	649.67	61.34	3	

**Delete selected peptide(s)**

Return to **Review proteins** locate the protein using the Search, as above, and clear the filter at the level of the peptides

There are now 7 remaining peptides for this protein

**Review selected protein**

Review the selected protein's identified peptides and validate their expression patterns.

**Accession:** gi|126700790  
**Description:** enolase (2-phosphoglycerate dehydratase) (2-phospho-D-glycerate hydro-lyase) [Clostridium difficile 630]

- Choose the level of detail**  
View the properties and expression profiles of either peptides or individual peptide ions:  
Show:  Peptides  Peptide ions  
Tip: you can also double-click a peptide to select and view its component ions.
- Compare expression profiles**  
Select peptides in the table to show their expression profiles in the chart below.
- Resolve any quantitative outliers**  
Tag any peptide ions whose expression profile is an outlier for this protein.  
[Learn about tagging and filtering](#)  
You can then review their identifications in more detail at the [Refine Identifications](#) step.

**Experiment design**  
Review your data from a different perspective:  
Current design: SAC

Identifier	Ions	Score	Correlation	Anova (p)	Max Fold Change	Highest Mean	Lowest Mean	Tag	Abundance	Neutral Mass
63.38_1327.7597n	1	88.370	0.971	0.00061	2.95	C	A		2.13E+05	1327.7597
50.10_1718.8842n	1	92.430	0.994	7.59E-05	2.95	C	A		2.842E+05	1718.8842
53.05_2188.0869n	2	106.390	0.983	0.00024	2.94	C	A		2.23E+05	2188.0869
39.09_1048.5911n	1	57.950	0.993	0.000319	2.78	C	A		1.052E+05	1048.5911
66.67_2046.0786n	2	105.600	0.900	0.00288	12.8	C	A		1.189E+05	2046.0786
41.67_1457.8020n	1	54.030	0.975	2.05E-05	5.12	C	A		1.996E+04	1457.802
37.60_1946.8902n	2	87.980	0.970	0.00118	4.49	C	A		3.744E+04	1946.8902

**Standardised Expression Profiles**

Standardised Normalised Abundance


Peptide profiles:  Contributes to protein measurements  Does not contribute to protein measurements

Now click **Return to the list of proteins** and create a tag for those proteins that have oxidised Methionine residues.


Modified proteins can be located by specifically searching for proteins containing modified peptides. Right click on the list of proteins and select **Modification** from the list of **Quick Tags**.

Accession	Peptides	Unique peptides	Confidence score	Anova (p)	q Value	Tag	Max fold change	Highest Mean	Lowest Mean
gi 296449145	4	4	450	0.086	0.00967		1.65	A	C
gi 255101963 (+2)	61	16	6.71E+03	0.000201	8.78E-05		2.52	A	C
gi 254976387 (+5)	63	18	7.07E+03	7.51E-05	4.33E-05		2.52	C	A
gi 1267004			81.7	0.000762	0.000203		1.72	C	A
gi 1267004			222	0.0108	0.00161		1.61	A	C
gi 2551019			641	---	---		---	---	---
gi 2551019				---	---		---	---	---
gi 2549740				---	---		---	---	---
gi 2549763				---	---		---	---	---
gi 254975791	2	2							
gi 254976383 (+2)	9	1							
gi 126699073	2	2	122	0.0143	0.00201		2.23	C	A
gi 126700407	16	16	1.64E+03	0.0498	0.00599		1.61	C	A

To find those proteins containing peptides with Oxidation on Methionine residues enter **Oxidation (M)**. This will automatically provide a named tag when you click **Create tag**.



To reduce the table to displaying only these proteins with modified peptides (Oxidation on methionine) use the tag filter to focus on these proteins.



The proteins table will now only display those proteins containing modified peptides.

QIP Progenesis QI.p Tutorial for DDA - Progenesis QI for proteomics

File Review Alignment Filtering Experiment Design Setup Review Peak Picking Peptide Ion Statistics Identify Peptides QC Metrics Refine Identifications Review Proteins Protein Statistics Report

nonlinear  
A Waters Company

Help

### Review Proteins

Using this screen, you can find the proteins of interest in your experiment.

**1 Set the quantitation options**  
If you've not already done so, choose between relative and absolute quantitation, use of Hi-N, protein grouping and more.  
[Protein options...](#)

**2 Create a shortlist to review**  
In the table, sort and filter the proteins based on their measurements, to generate a shortlist for further review.  
[How are the measurements calculated?](#)  
To sort the table by a given value, simply click the relevant column header.

**3 Review the proteins**  
For each protein of interest, review its peptide measurements and correlations:  
[View peptide measurements](#)  
You can also double-click to review a protein.

**4 Export data for further processing**  
By exporting your data to external tools, there's no limit to your analysis.  
[Export to pathways tool](#)  
[Export protein measurements](#)  
[Export peptide measurements](#)  
[Export peptide ion measurements](#)

**Experiment design**  
Review your data from a different perspective:  
Current design: AC

**Tag filter applied**  
proteins may be hidden [Edit...](#)

Accession	Peptides	Unique peptides	Confidence score	Anova (p)	q Value	Tag	Max fold change	Highest Mean	Lowest Mean	Description
gj 126699971	5	5	334	9.17E-06	1.47E-05	🌈	11.1	A	C	thioredoxin 2
gj 260682215 (+1)	24	11	2.08E+03	1.77E-05	2.1E-05	🌈	37.1	C	A	hemagglutinin,
gj 254973900 (+9)	15	15	1.95E+03	2.3E-05	2.51E-05	🌈	3.58	A	C	flagellin subun
gj 126697690	6	6	578	2.86E-05	2.8E-05	🌈	5.64	A	C	ferredoxin/flav
gj 126700790 (+1)	7	7	638	4.75E-05	3.3E-05	🌈	2.94	C	A	enolase (2-pho
gj 209571234	25	12	2.5E+03	4.78E-05	3.3E-05	🌈	5.99	A	C	cell wall protei
gj 254976387 (+5)	63	18	7.07E+03	7.51E-05	4.33E-05	🌈	2.52	C	A	cell surface pro
gj 126700129	3	3	267	0.000125	6.8E-05	🌈	2.32	A	C	translation inhi
gj 255101963 (+2)	61	16	6.71E+03	0.000201	8.78E-05	🌈	2.52	A	C	cell surface pro
gj 126697654	3	3	215	0.00023	9.14E-05	🌈	2.45	A	C	30S ribosomal
gj 126700372	2	2	126	0.00386	0.000743	🌈	4.95	A	C	PTS system HP
gj 54781347	6	6	544	0.00517	0.000903	🌈	1.9	A	C	2-hydroxyisoca
gj 126697631	7	7	626	0.00981	0.0015	🌈	1.41	A	C	50S ribosomal

**Selected protein: thioredoxin 2 (Trx2) [Clostridium difficile 630]**  
[View peptide measurements](#)

Quantifiable proteins displayed: 15

Section Complete




**Note:** the **Sequence Quick** tag can be used to locate Proteins containing peptides with specific sequence motifs, (i.e potential phosphorylation sites).

## Stage 13: Exporting Protein Data



Protein data can be exported in a csv file format. You can either export the **Protein, peptide or peptide ion measurements** using the options in the File Menu or use the buttons under Step 4 both available at the **Review Proteins** stage.

As an example of Data export use the Tag filtered set from the previous section for only proteins that have Oxidised Methionine residues.

First set the tag filter as shown below. Then select **Export Protein Measurements**.



The Export Protein Measurements dialog opens. Select the required fields and click OK. **Save** the file and then open the exported data file using the dialog that opens.




**Note:** if you have performed **Absolute Quantification** then the 'Amount' field will be available. Excel will open displaying the exported protein measurements

Accession	Peptide count	Unique peptides	q Value	Max fold change	Description	Normalized abundance					
						A1	A2	A3	C1	C2	C3
gi 126699971	5	5	1.33E-05	11.14094087	thioredoxin 2 (Trx2) [Clostridium difficile 6	1154.236	1215.294	1131.041	88.54863	109.4784	116.1807
gi 260682215;gi 209570719	24	11	1.78E-05	37.07319764	hemagglutinin/adhesin [Clostridium difficil	11.81321	16.35164	18.26058	530.6459	528.247	662.2461
gi 254973900;gi 260682017;gi	15	15	2.17E-05	3.581960969	flagellin subunit [Clostridium difficile QCD	5371.804	4705.512	5579.548	1446.399	1404.639	1519.994
gi 126697690	6	6	2.47E-05	5.642505943	ferredoxin/ferredoxin oxidoreductase subu	407.2672	466.0251	354.417	70.15977	74.2685	73.15401
gi 126700790;gi 296452046	7	7	3.00E-05	2.944231317	enolase (2-phosphoglycerate dehydratase)	62.67638	56.74192	64.69636	195.2659	177.3942	169.416
gi 209571234	25	12	3.00E-05	5.993442061	cell wall protein V [(Clostridium] difficile]	1917.599	1987.075	2256.879	369.1209	288.9127	370.0155
gi 254976387;gi 296452394;gi	63	18	3.86E-05	2.521368327	cell surface protein (S-layer precursor prote	25956.46	27703.39	25118.69	66993.88	60513.36	71122.46
gi 126700129	3	3	5.68E-05	2.318926178	translation inhibitor endoribonuclease [Clo	1468.906	1439.99	1571.524	672.1819	580.3695	679.5581
gi 255101963;gi 21702505;gi	61	16	7.81E-05	2.521286161	cell surface protein (S-layer precursor prote	69028.69	73480.05	61068.81	28145.33	24515.63	28082.57
gi 126697654	3	3	8.22E-05	2.451280888	30S ribosomal protein S8 [Clostridium diffic	256.4914	288.947	294.6438	119.5426	121.4714	101.6975
gi 126700372	2	2	0.000676	4.953785887	PTS system HPr protein [Clostridium difficil	26.49133	34.09065	24.56343	3.318793	7.936991	5.932161
gi 54781347	6	6	0.000887	1.901229315	2-hydroxyisocaproate-CoA transferase [Clo	399.6587	323.1495	309.6788	179.5084	155.1933	208.3611
gi 126697631	7	7	0.00149	1.412713841	50S ribosomal protein L7/L12 [Clostridium c	607.4534	520.8522	662.8525	440.9649	412.827	414.0928
gi 126700078	6	6	0.001572	1.398255682	molecular chaperone DnaK [Clostridium dif	130.7325	118.3408	103.9304	173.6002	165.4416	154.5476
gi 126697969;gi 255654423	10	10	0.0036	1.733236988	Beta-subunit of electron transfer flavoproti	1696.664	1586.885	2411.808	1149.335	897.5683	1239.062

**Note:** where there are multiple group members the other accession numbers are also exported.

At each stage in the Work flow there are a number of Export and Import options available from the **File** Menu. This includes the option to **Import Additional Protein Data** which can be used to increase the Protein meta data and also be used to sort the existing tabular data.



## Exporting Protein Data to Pathways Tool(s)

Using Progenesis, you can export protein lists to pathway analysis tools to help you understand your data in a wider biological context.

Currently Progenesis QI for proteomics supports the export to:

**IMPALA:** which aggregates and queries many other pathway analysis tools including KEGG and Reactome

**PANTHER:** classifies proteins for high-throughput analysis.


**MetaCore:** MetaCore is an integrated software suite for functional analysis of experimental data. MetaCore is based on a curated database of human protein-protein, protein-DNA interactions, transcription factors, signaling and metabolic pathways, disease and toxicity, and the effects of bioactive molecules.

**Ingenuity Pathways Analysis (IPA):** Explores experimental data by identifying relationships, mechanisms, functions, and pathways of relevance.

**Note:** Plugins for these tools are provided as standard.

**Note: Access to IPA** is dependent on having a local licence for IPA. Details of exporting data to IPA are described in Appendix 10 (page 93).

First select the protein data to export to the pathways tool using tag filtering to 'focus' the set to export. Then click **Export to pathways tool**, and select **IMPALA** as the tool followed by the test to be performed.



**Review Proteins**

Using this screen, you can find the proteins of interest in your experiment.

- Set the quantitation options**  
If you've not already done so, choose between relative and absolute quantitation, use of Hi-IL, protein grouping and more.
- Create a shortlist to review**  
In the table, sort and filter the proteins based on their measurements, to generate a shortlist for further review.
- Review the proteins**  
For each protein of interest, review its peptide measurements and correlations.
- Export data for further processing**  
By exporting your data to external tools, there's no limit to your analysis.

**Export Pathways Information**

Select a pathways tool  
Choose a pathways tool from the list below. You can find out more or download new plugins using the links below.

Which pathways tool do you want to use?

- IMPALA: Integrated Molecular Pathway Level Analysis
- IPA: Ingenuity Pathway Analysis
- MetaCore
- Panther Classification System

Accession	Max fold change	Highest Mean	Lowest Mean	Description
gi 123456789	182	C	A	cell surface proc...
gi 123456789	17.2	C	A	FMN-binding p...
gi 123456789	14.45	C	A	glyceraldehyde...
gi 123456789	13.4	C	A	ABC transporte...
gi 123456789	13.37	A	C	decarboxylase
gi 123456789	13.17	A	C	ribose-5-phosp...
gi 123456789	8.6	A	C	transketolase, c...
gi 123456789	8.54	C	A	fructose-1,6-bi...
gi 123456789	11.1	A	C	thioredoxin 2
gi 123456789	4.97	A	C	transcription el...
gi 123456789	4.98	C	A	cell surface proc...
gi 123456789	4.19	A	C	elongation fact...
gi 123456789	37.1	C	A	hemagglutinin...

Section Complete

Select either **Pathway over-representation analysis** or **Wilcoxon pathway enrichment analysis**.  
 Make sure the **Open IMPaLA in my browser** is ticked and then click **Copy proteins to clipboard**  
 When IMPaLA opens paste in the exported values and perform the test.

**Note:** currently the *Clostridium difficile* protein set being used in this user guide does not yield any pathway 'hits' with **IMPALA** or **Panther**.



Although the previous page outlines the process of exporting data to pathway tools, a text file (**Impala\_Wilcoxon\_demo.txt**) has been included in the tutorial download to allow the demonstration of exporting data to Impala for pathways analysis.

In Progenesis QI for proteomics set up the **Export to pathways tool** dialog to export the data to **IMPALA** to perform a Wilcoxon pathway enrichment analysis. Click copy to clipboard to open **IMPALA**.

To **simulate the Copy proteins to clipboard** open the text file in 'Notepad', select all and copy. In the **IMPALA** window, open in your browser, right click on the genes/protein panel and paste the contents of the file Impala\_Wilcoxon\_demo.txt.

Select UniProt as the identifier, tick Wilcoxon pathway enrichment analysis then Start Analysis.

**Note:** if you are using your own data then pressing **Copy proteins to clipboard** will open Impala and allow you to paste directly into Impala, without saving to a separate file.



Go back and start a new analysis

431 out of 2040 input gene identifiers were mapped to 431 distinct physical entities found in pathways.

**2008 pathways found.** Results per page: 50 of 41 (previous) (next)

download results

pathway name	pathway source	measured genes	all genes	P <sub>genes</sub>	Q <sub>genes</sub>
Gene Expression	Reactome	80	1251	0.000212	0.827
Metabolism of proteins	Reactome	86	693	0.0017	1
Alanine, aspartate and glutamate metabolism - Homo sapiens (human)	KEGG	13	35	0.00859	1
Ascorbate and aldarate metabolism - Homo sapiens (human)	KEGG	6	27	0.0313	1
Legionellosis - Homo sapiens (human)	KEGG	6	55	0.0355	1
glycolysis	HumanCyc	9	26	0.0391	1
Ribosome - Homo sapiens (human)	KEGG	46	135	0.0461	1
Metabolism of vitamins and cofactors	Reactome	8	88	0.0572	1
Glycerophospholipid biosynthesis	Reactome	7	96	0.0583	1
N-Glycan biosynthesis - Homo sapiens (human)	KEGG	5	49	0.0625	1
Histidine metabolism - Homo sapiens (human)	KEGG	5	23	0.0625	1
Fatty acid elongation - Homo sapiens (human)	KEGG	5	23	0.0625	1
PI3K-Akt signaling pathway - Homo sapiens (human)	KEGG	7	347	0.0625	1
Alzheimer's disease - Homo sapiens (human)	KEGG	17	168	0.0679	1
mRNA Splicing - Major Pathway	Reactome	12	131	0.0773	1
mRNA Splicing	Reactome	12	131	0.0773	1
Valine, leucine and isoleucine degradation - Homo sapiens (human)	KEGG	23	46	0.0774	1
RNA transport - Homo sapiens (human)	KEGG	10	165	0.084	1
Cholesterol biosynthesis	Reactome	4	22	0.0975	1

**Note:** in the case of exporting to **Panther** the normal process requires a file to be created and stored then opened in Panther to perform the pathway analysis.

Clear any applied tag filters then move to the next section/stage.




## Stage 14: Exporting identifications for submission to PRIDE


Following the analysis and review of your data you can export identification results as mzIdentML.

**Note:** this only supports identifications obtained from searches performed in Progenesis QI for proteomics version 4.0 or later, using the search methods **Mascot** or **Ion Accounting**.

Export Identifications is carried out at the Review Proteins stage of the workflow. Click on the **File** menu and select **Export mzIdentML for PRIDE submission...**




A dialog will appear allowing you to include the protein sequence or not in your export. Choosing to export the sequence will allow the display of the coverage map in PRIDE.



Enter a file name. Progenesis will save two files the mzIdentML as a .mzid file and also a corresponding .mgf file containing the spectra.


Open PRIDE Inspector and locate and select the mzid file to import.




Click open.

You will then be asked if you want to load the corresponding spectrum files. Click yes and locate the .mgf file.

A dialog opens and allows you to **Add spectra files**. Add the corresponding .mgf file and click **Set**



PRIDE Inspector will open allowing you to check the data.




Once checked you can submit the data to PRIDE using the PX Submission Tool.


**Note:** For exporting mzIdentML of results to PRIDE where the Search Method is Mascot, you **must select** additional **Mascot Export Fields** when exporting your search results from Mascot to Progenesis: (Protein sequence (Protein Hit Information) and Start and End (Peptide Match Information)), **Appendix 11 (page 96)**

## Stage 15: Protein Statistics

Protein Statistics opens with a Principal Components Analysis (PCA) for all the proteins displayed.




The Multivariate Stats can now be applied to all or subsets of proteins as determined by the current Tag filters (in this example none). Allowing you to identify similar patterns of expression using the Correlation Analysis.



Now move to the **Report** section to report on Proteins and /or peptides.


## Stage 16: Reporting

The **Report Design** stage allows you to select what views you want to include in a report based on the list of currently selected proteins.



As an example we will create a report for **only** the proteins containing peptides with Oxidation of Methionine residues.

1. First reduce the proteins to report on by selecting the '**Modification with Oxidation (M)**' tag. In this example it reduces the number of proteins in the table to 15.
2. Expand the various Report Design options (by default they are all selected)
3. Un-tick as shown below
4. Click **Create Report**



This opens a dialog to allow you to save the report, after which it will be opened in the form of a web page.

Click on the **Accession No.** in the proteins section of the Report and this will take you to the Assigned peptides for this protein

### Tutorial Methionine oxidised proteins

**Experiment:** Progenesis QI.p Tutorial for DDA  
**Report created:** 01/10/2019 10:08:29

**Proteins**

Protein building options  
 Protein grouping **Group similar proteins**  
 Protein quantitation **Relative Quantitation using Hi-3**

Accession	Peptides	Score	Anova (p)*	Fold	Tags	Description	Average Normalised Abundances	
							A	C
gi 254976387	63 (18)	7074.50	7.51e-005	2.52		cell surface protein (S-layer precursor protein) [Clostridium difficile QCD-66c26]	3.49e+007	8.80e+007
gi 255101963	61 (16)	6709.7						
gi 209571234	25 (12)	2502.3						
gi 260682215	24 (11)	2078.2						
gi 254973900	15 (15)	1945.7				flagellin subunit [Clostridium difficile QCD-66c26]		
gi 126697969	10 (10)	925.4						
gi 126700790	7 (7)	638.0						
gi 126697631	7 (7)	625.7						
gi 126700078	6 (6)	582.2						

**gi|254973900**

flagellin subunit [Clostridium difficile QCD-66c26]  
 15 peptides

Sequence	Peptide Ion	Score	Hits	Mass	Charge	Tags	Conflicts	Modifications	In quantitation	Average Normalised Abundances	
										A	C
AADDAAGLAISEK	148	100.52	5	1230.6087	2		0		yes	3.07e+006	1.00e+006
DTDVASEMVNLSK	1358	98.96	4	1407.6558	2		0		no	1.27e+005	3.41e+004
DTDVASEMVNLSK	3339	84.71	4	1423.6495	2		0	[8] Oxidation (M)	no	7.64e+004	1.48e+004
IADELLQLK											
IADELLQLKDEVER											
IADELLQLKDEVER											
IRDTDVASEMVNLSK											
IRDTDVASEMVNLSK											
IRDTDVASEMVNLSK											
IRDTDVASEMVNLSK											
ISSSIEFNGK											
LESTQNNLNNTIENVTAAESR											
LESTQNNLNNTIENVTAAESR											
MNILVQASQSMLAQANQQPQQ											

**Accession gi|254973900 (+9)**

**Description** flagellin subunit [Clostridium difficile QCD-66c26]  
**Peptides** 15 (15)  
**Score** 1945.78  
**Anova** 2.30e-005  
**Fold** 3.58

- Anova p-value ≤ 0.05
- Max fold change ≥ 2
- Modification with Oxidation (M)

**Note:** if you scroll down on the second page of the report you can locate expression profiles for each protein.


Having closed the report it can be reopened by double clicking on the saved html file.

**Note:** you can also copy and paste all or selected sections of the report to Excel and/or Word.

## Creating an Inclusion list

Inclusion lists can be used to try to increase the number of identified peptides you have in your experiment. They are used to control your mass spectrometer to try and concentrate the collection of MS/MS scans from specified m/z – retention time positions. Runs made using the inclusion list can then be imported into Progenesis QI for proteomics and the extra MS/MS scans added to the experiment.

As an example we are going to create an inclusion list for all the peptide ions that show a Significant difference between groups A and C (Anova  $p < 0.05$ ) and have **no MS/MS spectra**.



First return to **Review Peak Picking** using the Workflow icons.

#	Anova (p)	q Value	Fold	Tag	Notes	Highest Mean	Lowest Mean
10316	2.08E-13	4.84E-09	Infinity		Most abundant	C	A
9128	1.93E-12	2.24E-09	Infinity		Anova p-value ≤ 0.05	A	C
4026	3.08E-12	2.39E-09	Infinity		Max fold change ≥ 2	C	A
144	5.27E-12	3.07E-09	Infinity		Significantly up in C	C	A
7568	1.99E-11	9.25E-09	Infinity		Up regulated in C	A	C
6348	2.43E-11	9.42E-09	Infinity		Up regulated in A	A	C
1917	3.06E-11	1.02E-08	Infinity		Poorly Correlated	C	A
9494	5.34E-11	1.29E-08	Infinity			A	C
11023	6.02E-11	1.29E-08	Infinity				
314	6.35E-11	1.29E-08	Infinity				
3233	6.46E-11	1.29E-08	Infinity				
10463	6.94E-11	1.29E-08	Infinity				
1450	7.4E-11	1.29E-08	Infinity				

With no filters applied right click on a peptide ion in the table, select **Quick tags** and click on **No MS/MS data**.

**Create a filter**  
Show or hide peptide ions based on a selection of their tags. Move tags to the appropriate boxes to create the filter. For more guidance, please see the [online reference](#).

Available tags:

- Most abundant (2333 peptide ions)
- Anova p-value ≤ 0.05 (8174 peptide ions)
- Max fold change ≥ 2 (9822 peptide ions)
- Up regulated in A (4228 peptide ions)
- Significantly up in C (3946 peptide ions)
- Up regulated in C (3946 peptide ions)
- Poorly Correlated (2 peptide ions)

Show peptide ions that have all of these tags:

- No MS/MS data (8385 peptide ions)

Show peptide ions that have at least one of these tags:

Hide peptide ions that have any of these tags:

Clear the filter OK Cancel


Filter the table so that it is only showing peptide ions with **No MS/MS data**. To do this click **Create** and drag the tag on to the Show peptide ions that have all these tags and click OK.

The table will now only be displaying peptide ions that have no ms/ms.

#	Anova (p)	q Value	Fold	Tag	Notes	Highest Mean	Lowest Mean	m/z	z	Mass	RT (mins)	RT window (mins)	Abundance	Intensity	Max CV (%)	MS/MS
10316	2.08E-13	4.84E-10	Infinity			C	A	1002.0213	2	2002.028	33.698	0.212	8.82E+03	9.91E+04	0.73	0
9128	1.93E-12	2.24E-09	Infinity			A	C	592.2268	3	1773.659	42.473	0.294	5.93E+03	5.28E+04	1.22	0
144	5.27E-12	3.07E-09	Infinity			C	A	901.2218	2	1800.429	40.050	1.27	4.2E+06	1.43E+07	2.69	0
7568	1.99E-11	9.25E-09	Infinity		Add a note...	A	C	822.3762	2	1642.738	40.221	0.431	1.31E+04	7.11E+04	2.36	0
6348	2.43E-11	9.42E-09	Infinity			A	C	735.3522	3	2203.035	39.378	0.444	1.79E+04	9.02E+04	2.56	0
9494	5.34E-11	1.29E-08	Infinity			A	C	680.3218	3	2037.943	42.561	0.42	9.15E+03	4.99E+04	2.95	0
11023	6.02E-11	1.29E-08	Infinity			C	A	1446.7307	3	4337.170	53.154	0.227	2.2E+04	9.16E+04	3.27	0
3233	6.46E-11	1.29E-08	Infinity			C	A	533.9859	3	1598.936	31.074	0.455	1.87E+04	3.32E+05	3.28	0
10463	6.94E-11	1.29E-08	Infinity			A	C	1052.1714	3	3153.492	31.272	0.0772	8.33E+03	5.01E+04	3.1	0
4636	8.72E-11	1.35E-08	Infinity			A	C	666.6138	3	1996.820	43.337	0.633	3.15E+04	1.47E+05	3.7	0
5323	9.81E-11	1.43E-08	Infinity			A	C	946.7999	3	2837.378	34.444	0.48	4.09E+04	1.2E+05	3.94	0
5770	1.14E-10	1.46E-08	Infinity			A	C	1223.0393	2	2444.064	35.134	0.472	6.4E+04	1.07E+05	4.29	0
14169	1.16E-10	1.46E-08	Infinity			A	C	1012.2141	4	4044.827	29.887	0.264	3.67E+03	2.69E+04	3.22	0


Now select **Export Inclusion List...** from the file menu

Then select the appropriate format.



Finally export the file to an appropriate location

**Note:** with certain MS machines it is possible to widen the retention time windows being used, this can be controlled using the following dialog.



If you require further information on the inclusion list file formats then click the link **About this data format** in the Export Inclusion List dialog.

**Note:** The new LC-MS runs can then be added to the original experiment to increase the MS/MS coverage using the **Add files** facility at the Data Import Stage.

## Congratulations!

This document has taken you through a complete analysis using Progenesis QI for proteomics, from Alignment through the analysis workflow to generate lists of proteins exhibiting expression changes across biological conditions using powerful Multivariate Statistical analysis of the data.

Hopefully our example has shown you how this unique technology can deliver significant benefits with

- Speed
- Objectivity
- Statistical Power


If you would like to see the benefits of running Progenesis QI for proteomics using your own data and explore the Progenesis QI for proteomics workflow please go to Appendix 4: Licensing Runs (page 79).

## Appendix 1: Stage 1 Data Import and QC review of LC-MS data set

You can use your own data files, either by directly loading the raw files (Waters, AB SCIEX, Agilent, Bruker and Thermo) or, for other Vendors, convert them to mzXML or mzML format first.

To create a new experiment with your files: open Progenesis QI for proteomics and click **New**, bottom left of the **Experiments** page and give your experiment a name. Then select data type, the default is 'Profile data'.


**Note:** if you have converted or captured the data as centroided then select Centroided data and enter the Resolution for the MS machine used.



Click **Create experiment** to open the LC-MS Data Import stage of the workflow.

Select the 'Import Data file format', in this example they are mzXML files

Then locate your data files using **Import...**






Locate and select all the Data files (A1 to C3).

On loading the selected runs your data set will be automatically examined and the size of each file will be reduced by a 'data reduction routine', which reduces the data by several orders of magnitude but still retains all the relevant quantitation and positional information.

**Note:** For a large number of files this may take some time.



**Note:** you can start the automatic processing before the loading has completed.


For details of setting up the steps in the automatic processing wizard return to Stage 2A page 7.

## Additional details on Selection of Alignment Reference

### 1. Assess all runs in the experiment for suitability


- This method compares every run in your experiment to every other run for similarity.
- The run with the greatest similarity to all other runs is chosen as the alignment reference.

If you have no prior knowledge about which of your runs would make a good reference, then this choice will normally produce a good alignment reference for you. This method can take a long time.



### 2. Use the most suitable run from candidates that I select

- This method asks you to choose a selection of reference candidates, and the automatic algorithm chooses the best reference from these runs.



Select a subset of your runs as ‘alignment reference candidates’ and click **Mark selected runs**. An icon appears on the bottom left of each run to indicate that it is selected as a candidate.

This method is the most appropriate when you have some prior knowledge of your runs suitability as references:

- runs from pooled samples
- runs for one of your experimental conditions will contain the largest set of common peptides.

It is also the preferred option when analysing a large number of runs

### 3. Use this run

- This method allows you to manually choose the reference run.

Manual selection gives you full control, but there are a couple of risks to note:


- If you choose a pending run which subsequently fails to load, alignment will not be performed.
- If you choose a run before it fully loads, and it turns out to have chromatography issues, alignment will be negatively affected (for this reason we recommend that you let your reference run fully load and assess its chromatography before loading further runs).

### Review Chromatography

Each data file appears as a 2D representation of the run. If you created a **profile** experiment, at this stage you will be warned if any of the data files have been ‘centroided’ during the data acquisition and conversion process, as files must be of one format or the other.

You can delete run(s) by left clicking on the run in the list.

**Note:** you can also multi-select runs to remove by holding down the Ctrl key.



At the Import Data stage you can examine the quality of the imported runs using the 2D representation of the runs

**Import Data**

1 **Select your run data**  
Select one of the available data formats then click the Import button:  
Format:    
[About this data format](#) | [Download others](#)

2 **Perform automatic processing**  
Click the button below to start automatic processing. This will analyse your data as far as possible, before suggesting the next step.  
  
Automatic processing can be started while runs are importing.  
[Learn more about automatic processing](#)

Imported runs:

A1 A2 A3  
C1 C2 C3

**C1**

612.5 825.0 1037.5 1250.0 1462.5 1675.0 1887.5  
m/z

20  
40  
60  
80  
100  
Retention Time (min)

**About this run**

- MS peak count: 1,182,000
- MS/MS count: 9,542
- Total ion intensity: 3.643e+009
- Masked areas : none

**Note:** details of the current run appear on the top right of the view.

Once you have reviewed the imported runs click on Review Alignment or Section Complete to move forward to the Review Alignment Stage.


**Note:** you will be offered the automatic alignment if you have not performed the alignment automatically already.

Now move to the next stage in the workflow (page 14 in this user guide) by clicking **Section Complete**.


## Appendix 2: Stage 1 Processing failures

If a stage fails to complete successfully or only partially completes, the automatic processing dialog will warn you of the problem. This may or may not allow the automatic processing to complete.

For example: a run that fails to automatically align will trigger a warning, although analysis will continue; the automatic processing dialog will prompt you to 'drop-off' at the **Review Alignment** stage on completion to investigate the problem.




**Note:** in this example the run that failed to align will not contribute to the peak picking and will be excluded at the alignment stage (a cross appears in the include column).



You can either remove the run from the experiment at the **Import Data** or add it back in at the **Review Alignment** stage once the alignment of the run has been corrected.

As another example, runs that import successfully but with warnings at the **Import Data** stage will cause a flag in the readout to notify you of the potential quality issue.


If some runs fail to import (but not all), the automatic processing will continue, informing you that one or more runs have failed to import.



In this case you can remove the runs at Import Data and if appropriate replace them with additional runs.

**Note:** adding additional runs will then be aligned and peak picking should be re-done to include data from the added runs in the generation of the aggregate

An example of a problem that would halt the automatic processing would be the failure to successfully import all the potential reference candidates, (for example, where in importing, you specified selection of the alignment reference from several runs before they were fully imported and set the processing underway, and they later failed to import owing to problems with the runs).



In this case, the processing dialog would halt and prompt you to select another reference.


## Appendix 3: Stage 1 Data QC review and addition of exclusion areas

During the process of Data QC you may identify areas of the raw data for a particular run that appear 'noisy' yet still have identifiable 'isotopic patterns'.


For example if the run is part of a 'replicate set' of runs it is possible to exclude such areas on the noisy run by applying a mask to the area. By doing so this area is excluded during the initial part of the detection process in order that it does not 'interfere' with the detection of the peptide ions in the replicate group.

To do this select **Mask areas for peak picking** from the **Actions** menu on the top right of the screen.

Drag out an area over the noisy part of the run to create the mask.




**Note:** Click **Done** to return to the **Import Data** view, where you can zoom into the masked areas where you will see the isotopic peptide ions in the noise.



**Note:** if the level of noise is high and affecting many of your runs a preferred approach would be to re-optimize the chromatography to improve the levels of noise in your data.

## Appendix 4: Licensing runs (Stage 3)

When setting up a **New experiment** if you are evaluating Progenesis QI for proteomics with unlicensed runs then the licensing page will open after **Import Data** section.




**If you already have a programmed dongle attached to your machine then the License Runs page will not appear.**

To use this page to License your Runs you must first either obtain an **'Evaluation' Licence Code** from a **Sales Person** or **purchase a licence code directly**.

Each code will allow you to license a set number of runs.


The runs in your experiment will be listed as shown below.

To activate license(s) for the selected runs enter the code in the space provided and click **Use Licence code**.



**Note:** you will need an internet connection to use this method.

A message confirming successful installation of your licences will appear.



Click **OK**, the view will update and Alignment, the next stage in the workflow, will open with the licensed files.

## Appendix 5: Manual assistance of Alignment

### Approach to alignment

To place manual alignment vectors on a run (C2 in this example):

1. Click on Run C2 in the **Runs** panel, this will be highlighted in green and the reference run (in this example is C1) will be highlighted in magenta.
2. You will need approximately 5 **alignment vectors** evenly distributed from top to bottom of the whole run.
3. First drag out an area on the **Ion Intensity Map (C)**, this will reset the other 3 windows to display the same 'zoomed' area

**Review Alignment**  
Sample ions are aligned to compensate for drifts in retention time between runs.

1 **Align retention times automatically**  
For maximum reproducibility, the software can automatically align your runs.

2 **Review the alignment**  
Using the quality control measures, review and edit the runs' alignment:

1. Order the runs by alignment score and start by selecting the first run
2. Within each run, inspect and edit any areas rated as Needs Review

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/>	0	79.0%
A2	<input checked="" type="checkbox"/>	0	43.7%
A3	<input checked="" type="checkbox"/>	0	37.7%
C1	<input checked="" type="checkbox"/>	Ref	
C2	<input checked="" type="checkbox"/>	0	38.6%
C3	<input checked="" type="checkbox"/>	0	82.3%

Alignment quality: ■ Good ■ OK ■ Needs review

Section Complete


**Note:** the peptide ions moving back and forwards between the 2 runs in the **Transition** window (B) indicates the misalignment of the two runs.

**Note:** the Ion Intensity Map gives you a colour metric, visually scoring the current alignment. The overall score is placed next to the Vectors column in the table. With each additional vector added this score will update to reflect the 'changing' overall quality of the alignment. The colour coding on the Ion intensity Map will also update with each additional vector.


**Note:** The **Total Ion Chromatograms** window (D) also reflects the misalignment of the 2 runs for the current Retention Time range (vertical dimension of the current Focus grid in the **Ion Intensity Map** window).




- Click and hold on a green peptide ion in Window A as shown below



- As you are holding down the left mouse button (depending on the severity of the misalignment), the alignment vector will automatically find the correct lock. If not, drag the green peptide ion over the corresponding magenta peptide ion of the reference run. The red box will appear as shown below indicating that a positional lock has been found for the overlapping peptide ions.




**Tip:** while holding down the mouse button hold down the **Alt** key. This will allow smooth movement of the cursor as the **Alt** key allows you to override the 'automatic alignment' performed as you depress the mouse button.




- On releasing the left mouse button the view will 'bounce' back and a red vector, starting in the green peptide ion and finishing in the magenta peptide ion will appear.

**Note:** an incorrectly placed vector is removed by right clicking on it in the **Vector Editing** window


7. Now click **Show Aligned** on the top tool bar to see the effect of adding a single vector.




8. With the placement of a single manual vector the increase in the proportion of the **Ion Intensity Map** (C) showing green is reflected in the improved alignment score in the table. Now click in the Ion Intensity Map to relocate the focus in order to place the next manual vector.



9. Adding an additional vector will improve the alignment further as shown below.



10. The shift in the Retention Time (RT) is as a result of incorrect running of the chromatography. In many of these cases if the Automatic Alignment fails to generate a good alignment then removing all the alignment vectors for this run and placing a **single** manual vector to act as a 'seed' for the Automatic Alignment algorithm maybe all that is required to generate a good alignment.
11. In the case of the example shown above placing 2 vectors from the top to the bottom of the run is sufficient to markedly improve the alignment (**Note:** the improvement in the score with the addition of only 2 vectors).
12. At this point you would redo the automatic alignment of this image by selecting **Align runs automatically**. **Note:** if you are focusing only on the alignment of one run, then un-tick the other runs in the alignment dialog.



The number of manual vectors that you add at this stage is dependent on the misalignment between the current run and the Reference run.

**Note:** In many cases only using the Automatic vector wizard will achieve the alignment.

### Review Alignment

Sample ions are aligned to compensate for drifts in retention time between runs.

[Learn about the visualisations shown here](#)

- Align retention times automatically**  
For maximum reproducibility, the software can automatically align your runs.
- Review the alignment**  
Using the quality control measures, review and edit the runs' alignment:
  - Order the runs by alignment score and start by selecting the first run
  - Within each run, inspect and edit any areas rated as Needs review

[Learn about the review and editing process](#)

Run	Include?	Vectors	Score
A1	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	0	79.0%
A2	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	0	43.7%
A3	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	0	37.7%
C1	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	Ref	
C2	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	355	97.9%
C3	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/>	0	82.3%

Ion maps: ■ Alignment target ■ Run being aligned

Show Aligned
Show Unaligned
Remove Vectors

#### Vector editing

#### Transition

#### Ion intensity map

#### Total ion chromatogram

Alignment quality: ■ Good ■ OK ■ Needs review

**Tip:** a normal alignment strategy would be: to run the automatic alignment first for all runs, then order the alignments based on score. For low scoring alignments remove all the vectors and place 1 to 5 manual vectors to increase the score then perform automatic alignment. Then review the improved alignment score.

To review the vectors automatic and manual, return to page 15.

## Appendix 6: Within-subject Design

To create a **Within-subject Design** for your data set select this option on the **Experiment Design Setup** page and enter the name of the design.

In this example there are 3 Subjects (i.e. patients A, B and C) who have been individually sampled: Before(1), During (2) and After (3) treatment

Which experiment design type do you want to use for this experiment?

Between-subject Design

Do samples from a given subject appear in only one condition? Then use the between-subject design.

To set up this design, you simply group the runs according to the condition (factor level) of the samples. The ANOVA calculation assumes that the conditions are independent and therefore gives a statistical test of whether the means of the conditions are all equal.

**Within-subject Design**

Have you taken samples from a given subject under different conditions? Then use the within-subject design.


A standard ANOVA is not appropriate because the data violates the ANOVA assumption of independence. With a repeated measures ANOVA individual differences can be eliminated or reduced as a source of between condition differences (which helps to create a more powerful test).

The within-subject design can be thought of as an extension of the paired-samples t-test to include comparison between more than two repeated measures.

	Before	During	After
Patient X	X1	X2	X3
Patient Y	Y1	Y2	Y3
Patient Z	Z1	Z2	Z3

When the design page opens use the **Add Subject** and **Add Condition** buttons to create the matrix that fits your experimental design, over typing the names as required.

Then click on 'Select Sample' in each box of the matrix and select the appropriate sample.






You can create additional Experimental Designs using the New tab

All of these Experimental Designs are available at the later stages of the workflow with the exception of **Identify Peptides** and **Refine Identifications** (including Resolve Conflicts).

## Appendix 7: Power Analysis (Peptide Ion Stats)

Power analysis is a statistical technique that is used to gauge how many replicates are needed to reliably see expression differences in your data. It is available through the Peptide Ion Statistics section of the workflow.

To perform a power analysis of the data click on **Ask another question** at the top of the table in the Peptide Ion Statistics section. A selection of 3 tools will appear.

	<b>Principal Components Analysis</b> Are there any outliers in my data? Does my data cluster according to my experimental conditions?
	<b>Correlation Analysis</b> Group my peptide ions according to how similar their expression profiles are.
	<b>Power Analysis</b> How many replicates should I run? What is the power of my experiment?


Select the option

*'How many replicates should I run and what is the power of my experiment?'*

It answers this question by informing you:

*'How many replicates you need so that at least 80% of your peptide ions with a power >0.8'*

Using the **Anova p-value  $\leq 0.05$  peptide ions (8174)**, as an example, view the power analysis.



This is displayed graphically showing that 79.1% of the 8174 peptide ions have a power of 80% or that 4 replicates would give you 96.7% of your data with power > 0.8.

- **The power of a statistical test** reflects our confidence in the experimental data's ability to find the differences that do actually exist
- The power is expressed as a percentage, where 80% power is an accepted level, therefore allowing you to assess the number of sample replicates that would be required to achieve a power of 80%.

## Appendix 8: Resolve Conflicts

This stage allows you examine the behaviour of the identified peptides and resolve any conflicts for the various peptide assignments at the protein level.

The **Resolve Conflicts** stage (now accessed at the bottom left of the Refine Identifications stage) provides a number of interrelated graphical and tabular views to assist you in the manual validation of the peptides that have been assigned to proteins and also to review the relevance of the data returned from the search.

Progenesis QI for proteomics

File Import Data Review Alignment Filtering Experiment Design Setup Review Peak Picking Peptide Ion Statistics Identify Peptides QC Metrics Refine Identifications Review Proteins Protein Statistics Report

Experiment design: AC

Proteins No filter applied Create...

Peptide ions of gi|254976387 No filter applied Create...

Accession	Peptides	Unique	Conflicts	Score
gi 254976387 (+5)	63	18	88	7.07E+03
gi 255101963 (+2)	61	16	88	6.71E+03
gi 209571234	25	12	18	2.5E+03
gi 260682215 (+1)	24	12	18	2.08E+03
gi 126700407	16	16	0	1.64E+03
gi 254973900 (+9)	15	15	0	1.95E+03

#	Σ	Score	Hits	Correlation	Mass	Mass error (p...)	RT (mins)	Charge	Tag	Abundance	Confli
1	2	75.2	5	-0.624	2413.301	2.48	55.6	3	Red	1.24E+08	1
17	0	128	5	0.464	2413.296	0.265	55.6	2	Red	9.38E+07	1
533	0	63.7	2	-0.931	2413.295	-0.243	55.6	4	Blue	1.25E+06	1
761	0	61.3	4	-0.791	2413.295	0.0288	50.5	3	Red	3E+05	1
1256	2	123	5	0.844	2413.294	-0.699	50.5	2	Blue	2.84E+05	1
2421	0	122	5	0.494	2413.294	-0.603	53.1	2	Red	1.67E+05	1

Protein: gi|254976387 cell surface protein (S-layer precursor protein) [Clostridium difficile QCD-66c26]

Peptide ions of selected protein

Peptide ion views Protein resolution

Standardised Expression Profiles

Peptide ions of selected protein

Refine Identifications Protein options... Section Complete

The number of conflicts you have to resolve will depend on the scope and stringency of the filters you apply at the **Refine Identifications** stage.

**Note:** the default **Protein options** for protein grouping and Protein quantitation are set as shown

QI Protein quantitation options

Quantitation method:  
Relative Quantitation using Hi-N

Number of peptides to measure per protein (N): 3

Employ protein grouping, i.e. hide proteins whose peptides are a subset of another protein's.

OK Cancel

**Note:** manual conflict resolution is not essential to performing quantitation however it can lead to more quantifiable proteins.


This means that if you choose **not** to resolve the conflicts then proteins, to be considered for quantitation, require at least one unique peptide (number in brackets after peptide count).

For more details on **Protein Grouping** go to page 49 and **Protein Quantitation** go to page 51


Open **Resolve Conflicts** and order the data in the Proteins table (A) on the basis of **Conflicts**.

**Note: the look of the tables (with regards to ordering) in the following section may vary slightly.**

Make 'cell surface protein (S-layer precursor etc ' the current protein by clicking on it in Window A (a circular orange symbol indicates current protein). This protein has 63 peptides assigned (window B) which have a total of 88 conflicts. To view the conflicting assignments click on the **Protein Resolution** tab (window C previous page).



In this case the conflicting peptide assignments are with 'The same protein' (from a different strain) which also contains 88 conflicts. A simple resolution to these conflicts is to right click on the conflicting protein and turn off all its peptides (based on lower number of peptides and score).



**Note: as you un-assign the peptides the number of conflicts update 'on the fly' in all the windows.**



In this case the conflicting peptides are unassigned from the 'precursor' protein.

In some cases you can resolve the conflicts between 2 proteins on the basis of consistent peptide expression. In the example below the proteins share 18 conflicts. For the protein showing 2 clear patterns of expression you can un-tick all the peptides with conflicts in the corresponding peptides table

Progenesis QI.p Tutorial for DDA - Progenesis QI for proteomics

Experiment design: AC

Proteins: No filter applied

Accession	Peptides	Unique	Conflicts	Score
gi 254976387 (+5)	63	63	0	7.07E+03
gi 255101963 (+2)	0	0	0	0
gi 209571234	25	12	18	2.5E+03
gi 260682215 (+1)	24	11	18	2.08E+03
gi 126700407	16	16	0	1.64E+03
gi 254973900 (+9)	15	15	0	1.95E+03

Peptide ions of gi|209571234

#	Σ	Score	Hits	Correlation	Mass	Mass error (p...)	RT (mins)	Charge	Tag	Abundance	Conflict
2723	0	89.9	5	-0.461	2408.236	-0.344	51	2	1.04E+05	1	
206	1	98.9	5	-0.607	1502.794	-0.201	29.7	2	8.33E+05	1	
2974	0	---	---	0.063	1502.793	-0.773	29.7	3	1.43E+04	1	
234	1	41.7	5	-0.575	1732.899	-0.448	30	3	6.87E+05	1	
484	0	77.3	5	-0.596	1732.899	-0.528	30	2	5.26E+05	1	
284	0	52.3	5	-0.542	1051.555	-0.192	24.6	2	4.97E+05	1	

Protein: gi|209571234 cell wall protein V [[Clostridium] difficile]

Peptide ions of selected protein

Standardised Expression Profiles

Peptide ions of selected protein

This leaves the peptides with the same expression pattern assigned with the appropriate protein thus resolving the conflicts.

Progenesis QI.p Tutorial for DDA - Progenesis QI for proteomics

Experiment design: AC

Proteins: No filter applied

Accession	Peptides	Unique	Conflicts	Score
gi 254976387 (+5)	63	63	0	7.07E+03
gi 255101963 (+2)	0	0	0	0
gi 209571234	12	12	0	1.25E+03
gi 260682215 (+1)	24	24	0	2.08E+03
gi 126700407	16	16	0	1.64E+03
gi 254973900 (+9)	15	15	0	1.95E+03

Peptide ions of gi|209571234

#	Σ	Score	Hits	Correlation	Mass	Mass error (p...)	RT (mins)	Charge	Tag	Abundance	Conflict
1526	1	67	5	---	1609.903	-0.277	61.6	3	5.89E+04	0	
910	1	61.4	4	---	1115.571	0.124	28.5	2	1.44E+05	0	
2165	1	108	4	---	1821.024	-0.65	38.7	2	9.46E+04	0	
437	0	93.7	4	1.000	1588.82	-0.343	40.2	2	8.44E+05	0	
411	1	116	5	1.000	1708.671	2.41	46.9	2	1.52E+06	0	
727	1	48.3	5	0.998	1836.96	-1.35	39.1	3	4.29E+05	0	

Protein: gi|209571234 cell wall protein V [[Clostridium] difficile]

Peptide ions of selected protein

Standardised Expression Profiles

Peptide ions of selected protein

**Note:** the number of conflicts you have to resolve will depend on the scope and stringency of the filters you apply at the **Refine Identifications** stage.

To resolve the remaining conflicts in the example shown below first order the conflicts in panel B and select the top one (which may still be selected) Panel B will display the peptides for this protein and the number of conflicts for each peptide. Panel E will also update to show the conflicting protein.

The screenshot shows the Progenesis QI interface with several panels. The top panel displays a list of proteins with columns for Accession, Peptides, Unique, Conflicts, and Score. The second panel shows a detailed view of peptide ions for a selected protein, with columns for #, Score, Hits, Correlation, Mass, Mass error, RT, Charge, Tag, Abundance, and Conflict. Annotations include:

- B**: A box around the peptide ion entry for mass 1198.656.
- E**: A box around the protein entry for *gj|254976383* in the 'Conflicting proteins for peptide ion 3657' table.
- F**: A box around the peptide ion entry for mass 2319.16 in the 'Peptide ions of *gj|10180205*' table.

Favouring the protein with the higher score, resolve the conflict by switching off (or un-assigning) the peptide in panel F for the protein with the lower score. By doing this, the other 3 panels update to show the change in conflicts.

This screenshot shows the interface after resolving a conflict. The 'Conflicts' column in the protein table is updated. In the 'Peptide ions of *gj|10180205*' table, the entry for mass 3098 now has a conflict count of 0. The 'Conflicting proteins for peptide ion 3098' table shows that *gj|10180205* is now the only protein associated with this peptide.

Continue with **Conflict resolution** until no remaining conflicts in the Proteins table.

This screenshot shows the next step in conflict resolution. The 'Conflicts' column in the protein table is updated again. The 'Conflicting proteins for peptide ion 165' table shows that *gj|126699078* is now the only protein associated with this peptide.

**Note:** the abundances will need to be recalculated as a result of performing Conflict resolution. This is achieved by clicking on the Recalculate abundances, which appears during Conflict resolution

A close-up view of the bottom right corner of the interface, showing the 'Recalculate abundances' button highlighted with a red box.


## Appendix 9: Using Clip Gallery to Save and Export Pictures and Data

At every stage of the Progenesis QI for proteomics workflow the images and data tables can be added to the Clip Gallery.

The saved images are retained as part of the experiment and are stored accordingly. This facility allows you to capture (high resolution) images that can be used in the development of specific reports and/or used as part of the process of publishing your experimental findings.


As an example of using the Clip Gallery, at the **Protein Statistics** view, displaying the PCA plot.

Right click on the **Biplot** View and select Add to clip gallery....




This will open a dialog displaying what is to be saved and allows you to alter the title and provides a description of the image for later reference.


Enter details as required and click **Add to clip gallery**



To view, edit and/or export from the clip galley the gallery can be accessed from the **File** menu.



Selecting an image in the gallery makes available an **Actions** menu that allows you to manipulate the output of the image.



**Note:** there is also the capacity to **Export all...** the images in the experiments clip gallery which creates a list of files in a folder of your choice where the file name is based on the image title.

**Note:** right clicking on a table and adding it to the Clip Gallery allows you to export the current content to Excel.

## Appendix 10: Exporting protein identities to IPA

**IPA (Ingenuity Pathway Analysis)**: is a commercial web-based software application for the analysis, integration and interpretation of data derived from any experiments that generate gene and chemical lists with expression data, such as metabolomics and proteomics.

Before exporting data to IPA, make sure that you can launch IPA using the default browser on your system. The plug-in relies on the browser to upload the Progenesis-exported data to IPA.

**Note:** the use of IPA requires a licensed copy, for licensing details please visit the [IPA website](#).

First select the Protein data to export to the pathways tool using tag filtering to 'focus' the set to export.

Click **Export to pathways tool** and select IPA and click **Next**.

The screenshot shows the Progenesis QI software interface. The main window is titled "QIP Progenesis QI.p Tutorial for DDA - Progenesis QI for proteomics". The "Review Proteins" step is active, showing a table of protein data with columns: Accession, Peptides, Unique peptides, Confidence score, Anova (p), Tag, Max fold change, Highest Mean, and Lowest Mean. A dialog box titled "QIP Export Pathways Information" is open, asking "Which pathways tool do you want to use?". The dropdown menu shows "IPA: Ingenuity Pathway Analysis" selected. Below the dialog, there are buttons for "< Back", "Next >", and "Cancel".

Accession	Peptides	Unique peptides	Confidence score	Anova (p)	Tag	Max fold change	Highest Mean	Lowest Mean
gi 54781345 (+1)	5	5	406	1.13E-05		11.1	A	C
gi 54781347	7	7	622	0.0042		1.93	A	C
gi 126697687	2	2	137	0.000376		3.38	A	C

Select the appropriate identifier type used for proteins in the experiment.

**Note:** for proteins with different identifier types, filter the protein list to contain only proteins with single identifier type at a time, then perform the export operation and merge the resulting lists in IPA.

The screenshot shows the "QIP Export Pathways Information" dialog box. The "Configure your export" section is active, asking "Choose which identifications to export and the type of analysis you want to perform." The dropdown menu for "Select the identifier type used for proteins in the experiment:" shows "UniProt/Swiss-Prot Accession" selected. Below the dropdown, there is a yellow warning box: "For an overrepresentation analysis, filter your data to only show proteins with significantly altered expression." At the bottom, there are buttons for "< Back", "Export proteins to IPA", and "Close".

Finally select the way you want to export the proteins to IPA.

If you are performing over-representation analysis or other types of analysis that do not require the protein expression data then you can make use of the **Create a list** option to export data to IPA.

**Note:** using this option you will be asked to filter your data to only show proteins with significantly altered expression

If you require to perform enrichment analysis or other types of analyses that use protein expression data, select the option to **Upload expression dataset** to IPA.

**Note:** You need to select the type of protein identifier used by at least one protein whose expression data you want to export. If there are proteins with different identifier types, you will be able to select additional identifier types in IPA after the dataset is uploaded.

The operation exports data from a single observation at a time - select the two experimental conditions (taken from the current **experiment design as set in Progenesis QI for proteomics** you wish to compare as that observation).

Click **Export proteins to IPA** .

If you have access to a licenced copy of IPA then log in

IPA will open displaying the imported data from Progenesis QI for Proteomics.

The screenshot displays the IPK IPA software interface. At the top, there is a menu bar (File, Edit, View, Window, Help) and a search bar for genes and chemicals. Below this, the 'Dataset Upload - New Dataset 2016-05-20 10:29 AM' section is active, showing configuration options for file format (Flexible Format), column headers, identifier type (UniProt/Swiss-Prot Accession), and array platform (Not specified/applicable). A 'Raw Data (512)' table is visible, listing protein IDs and their corresponding expression values. To the right, the 'Dataset Upload Workflow Instructions' panel provides a step-by-step guide for uploading data, including a 'Data Upload Workflow' section and a 'GenBank' dropdown menu.

ID	Observation 1	Ignore	Ignore	Ignore	Ignore	Ignore
1	P34731	-1.00902155389087	---	---	---	---
2	C4YQR7	1.03618874042727	---	---	---	---
3	C4YMC3	-1.1630436093987	---	---	---	---
4	C4YQN7	-1.14248088552437	---	---	---	---
5	P43098	-1.00167260838215	---	---	---	---
6	C4YR46	1.11144195146305	---	---	---	---
7	O13430	1.0005161363124	---	---	---	---
8	P46587	1.16697510036399	---	---	---	---
9	P82610	1.14143675382931	---	---	---	---
10	P46598	1.06004923965961	---	---	---	---
11	Q96V89	-1.10613636924535	---	---	---	---
12	P41797	1.50747573107542	---	---	---	---
13	P28877	-1.20624862059398	---	---	---	---
14	C4YK39	1.05494661189734	---	---	---	---
15	Q59KZ1	1.03999813845519	---	---	---	---
16	C4YL05	-1.0312270752554	---	---	---	---
17	P46273	1.31165982609058	---	---	---	---
18	C4VIL8	1.26909493012736	---	---	---	---
19	O94039	1.55325257659893	---	---	---	---
20	O13287	1.38535524141575	---	---	---	---

You can now explore your protein expression data using the tools available in IPA. Instructions on how to manage and explore your imported data set are provided by IPA.

## Appendix 11: Exporting protein identities to MetaCore

**Metacore:** is a commercial web-based software application for the analysis, integration and interpretation of data derived from any experiments that generate gene and chemical lists with expression data, such as metabolomics and proteomics.

Before exporting data to MetaCore, make sure that you can launch MetaCore using the default browser on your system. The plug-in relies on the browser to upload the Progenesis-exported data to MetaCore.

First select the Protein data to export to the pathways tool using tag filtering to 'focus' the set to export.

Click **Export to pathways tool** and select MetaCore

The screenshot shows the Progenesis QI software interface. The main window is titled 'QIP Progenesis QI Tutorial for DDA - Progenesis QI for proteomics'. The 'Review Proteins' section is active, showing a table of protein data with columns for Accession, Peptides, Unique peptides, Confidence score, Anova (p), q Value, Tag, Max fold change, Highest Mean, Lowest Mean, and Description. A dialog box titled 'QIP Export Pathways Information' is open, asking 'Which pathways tool do you want to use?'. The 'MetaCore' option is selected in the dropdown menu. The 'Export to pathways tool' button is highlighted in red. The 'Section Complete' button is visible at the bottom right.


Then click **Next** to Configure your export page .

Select the appropriate identifier type used for proteins in the experiment.

The screenshot shows the 'QIP Export Pathways Information' dialog box. The 'Configure your export' section is active, asking 'Choose which identifications to export and the type of analysis you want to perform.' The 'Select the identifier type used for proteins in the experiment:' dropdown menu is open, showing options: Swiss-Prot Accession, GI Number, ENSEMBL, and EMBL. 'Swiss-Prot Accession' is selected. The 'Export proteins to MetaCore' button is highlighted in blue.




If you have multiple conditions in your current experiment design you can select which condition should be used as the control



Finally click **Export proteins to MetaCore**.

If you have not already done so, you will be prompted to enter your MetaCore *User name* and *Password*; once you have done this, your data will be exported to your MetaCore account, and the MetaCore user interface will open in your default web-browser.



Metacore will open displaying the imported data from Progenesis QI for Proteomics.

You can now explore your protein expression data using the tools available in MetaCore. Instructions on how to manage and explore your imported data set are provided by MetaCore.

## Appendix 12A Search engine parameters for Mascot

The parameters applied to the Mascot search that yielded the search results used in this user guide are shown below:

### MASCOT MS/MS Ions Search

<b>Your name</b> <input type="text" value="andy b"/>	<b>Email</b> <input type="text" value="andy.borthwick@nonlinear.com"/>
<b>Search title</b> <input type="text" value="Progenesis QI search"/>	
<b>Database(s)</b> <div style="border: 1px solid gray; padding: 2px;">           Mark_Test2            MSDB  <b>NCBIInr</b>            NCBIInr_UserGuide            NIST_Yeast_Sigma         </div>	<b>Enzyme</b> <input type="text" value="Trypsin"/>
<b>Taxonomy</b> <input type="text" value="... Firmicutes (gram-positive bacteria)"/>	<b>Allow up to</b> <input type="text" value="1"/> missed cleavages
<b>Fixed modifications</b> <div style="border: 1px solid gray; padding: 2px;">           Carbamidomethyl (C)         </div> <input type="checkbox"/> Display all modifications	<b>Quantitation</b> <input type="text" value="None"/>
<b>Variable modifications</b> <div style="border: 1px solid gray; padding: 2px;">           Oxidation (M)         </div>	<div style="border: 1px solid gray; padding: 2px;">           Acetyl (K)            Acetyl (N-term)            Acetyl (Protein N-term)            Amidated (C-term)            Amidated (Protein C-term)            Ammonia-loss (N-term C)            Biotin (K)            Biotin (N-term)            Carbamyl (K)            Carbamyl (N-term)            Carboxymethyl (C)         </div>
<b>Peptide tol. ±</b> <input type="text" value="9"/> ppm <input type="text" value="# 13C 0"/>	<b>MS/MS tol. ±</b> <input type="text" value="0.6"/> Da
<b>Peptide charge</b> <input type="text" value="2+"/>	<b>Monoisotopic</b> <input checked="" type="radio"/> Average <input type="radio"/>
<b>Data file</b> <input type="text" value="C:\Users\andy.borthwick\Document"/> <input type="button" value="Browse..."/>	<b>Precursor</b> <input type="text"/> m/z
<b>Data format</b> <input type="text" value="Mascot generic"/>	<b>Error tolerant</b> <input type="checkbox"/>
<b>Instrument</b> <input type="text" value="ESI-TRAP"/>	<b>Report top</b> <input type="text" value="AUTO"/> hits
<b>Decoy</b> <input type="checkbox"/>	
<input type="button" value="Start Search ..."/>	<input type="button" value="Reset Form"/>

Database : NCBIInr (circa 04/16) was used with the Taxonomy restriction set to Firmicutes

Fixed modifications: Carbamidomethyl (C) and variable modification Oxidation (M)

Peptide Tol: 9ppm

Instrument: ESI-Trap

## Appendix 12B Use Additional Export Fields in Mascot for PRIDE

For exporting mzIdentML of results to PRIDE from Progenesis QI for proteomics you **must select** additional Mascot Export Fields when exporting your search results from Mascot: (Protein sequence (Protein Hit Information) and Start and End (Peptide Match Information))

<b>Protein Hit Information</b>	<input checked="" type="checkbox"/>
Score	<input checked="" type="checkbox"/>
Description *	<input checked="" type="checkbox"/>
Mass (Da) *	<input checked="" type="checkbox"/>
Number of queries matched	<input checked="" type="checkbox"/>
Percent coverage **	<input type="checkbox"/>
Length in residues **	<input type="checkbox"/>
pI **	<input type="checkbox"/>
Taxonomy **	<input type="checkbox"/>
Taxonomy ID **	<input type="checkbox"/>
Protein sequence **	<input checked="" type="checkbox"/>
emPAI	<input type="checkbox"/>
<small>* Occasionally requires information to be retrieved from external utilities, which can be slow</small>	
<small>** Always requires information to be retrieved from external utilities, which can be slow</small>	
<b>Peptide Match Information</b>	<input checked="" type="checkbox"/>
Experimental Mr (Da)	<input checked="" type="checkbox"/>
Experimental charge	<input checked="" type="checkbox"/>
Calculated Mr (Da)	<input checked="" type="checkbox"/>
Mass error (Da)	<input checked="" type="checkbox"/>
Start	<input checked="" type="checkbox"/>
End	<input checked="" type="checkbox"/>
Number of missed cleavages	<input checked="" type="checkbox"/>

## References

Silva, JC, Gorenstein, MV, Li, G-Z, Vissers, JPC and Geromanos. Molecular and Cellular Proteomics (2006); 5 : 144-156 Absolute Quantification of Proteins by LCMS<sup>E</sup>